



Cytek® Northern Lights™ -CLC User's Guide



Copyrights

© 2020-2023, Cytek Biosciences Inc. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without prior written permission from Cytek Biosciences.

The information in this guide is subject to change without notice. Cytek Biosciences reserves the right to change its products and services at any time to incorporate the latest technological developments. Although this guide has been prepared with every precaution to ensure accuracy, Cytek Biosciences assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information. Cytek Biosciences welcomes customer input on corrections and suggestions for improvement.

Trademarks

Cytek, the Cytek logo, SpectroFlo, Northern Lights, cFluor, and Complexity and Similarity are trademarks of Cytek Biosciences, Inc. All other service marks, trademarks and tradenames appearing in this user guide are the property of their respective owners.

FCC Information

WARNING: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTICE: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, can cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his or her own expense.

Shielded cables must be used with this unit to ensure compliance with the Class A FCC limits.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

Cet appareil numérique de la classe A respecte toutes les exigences du Règlement sur le matériel brouilleur du Canada.

CDRH Information

Class I laser product.

Regulatory Information

The Northern Lights™-CLC flow cytometer with SpectroFlo® CLC software application is CE marked in compliance with Regulation EU 2017/746 on In vitro Diagnostic Medical Devices. This device is intended for use with in vitro diagnostic flow cytometric applications, only in countries where the regulation approval has been obtained from the local regulatory authorities.

Class 1 Laser Product



For In Vitro Diagnostic Use

Authorized Representative in the European Community



Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

History

Revision	Date	Change
N9-20039US Rev. A	7/2020	Initial release
N9-20039US Rev. B	1/2021	Added information on new loader
N9-20039US Rev. C	9/2021	Minor corrections per reviewers
N9-20039US Rev. D	1/2022	Corrected addresses and added manufacture symbol. Updated Replacement part numbers table to align with current released parts and assemblies.
N9-20039US Rev. E	11/2023	Updates to replacement parts and specifications tables.

Contents

Chapter 1: Introduction	9
Intended Purpose	9
About this Guide	9
Safety	9
Safety Symbols	9
General Safety	11
Electrical Safety	12
Biological Safety	12
Laser Safety	12
Technical Support	13
Chapter 2: Overview	15
Cytometer Overview	16
Fluidics System	17
Optics	19
Loader Overview	20
Software Overview	22
Get Started Menu	22
About Experiments	23
About Worksheets	26
About Reports	28
Chapter 3: Startup & Shutdown	29
Filling the Sheath and Emptying the Waste	29
Filling the Sheath	29
Emptying the Waste	30
Starting Up the System	31
Shutting Down the System	32
Chapter 4: QC & Setup	33
Daily QC	33
Performing Daily QC	33
QC Report	36
Reference Controls	38
Running Reference Controls	38
Updating Reference Controls	43
Levey-Jennings Tracking	44
Gain Settings	45

Alarm Ranges	45
Chapter 5: Acquisition	47
Raw vs Unmixed Data	47
Unmixing and Compensation	47
Setting Up an Experiment	48
Acquisition Experiment Overview	48
Experiment Display	48
Creating a Default Experiment	57
Creating a New Experiment	58
Creating an Assay	64
Create New Assays	64
Chapter 6: Library, Preferences, and Users	67
Library	67
QC Beads	67
Fluorescent Tags	67
Labels	68
Keywords	69
Spillovers	71
User Settings	71
Loader Settings	72
Worksheet Templates	72
Experiment Templates	73
Backup & Restore	73
Preferences	75
Acquisition Preferences	75
Worksheet Preferences	77
Plot Preferences	78
Gates Preferences	79
Statistics Preferences	80
Fonts Preferences	81
Notifications Preferences	82
Storage Preferences	82
QC Setup Preferences	84
Cytometer Preferences	84
Users	86
Managing Users	86
Use Time	88
User Role	88
User Policy	91
Chapter 7: Unmixing and Compensation	93
Spectral Unmixing	93
Understanding Full Spectrum Flow Cytometry	93

Unmixing Workflows	94
Unmixing Overview	94
Live Unmixing	96
Post-Acquisition Unmixing	101
Virtual Filters	106
Chapter 8: Worklist	111
Create or open a worklist	111
Add tasks	112
Reference controls	113
Unmix	114
Lot specific reference controls	115
Arrange sample locations in a carrier	116
Data collection	116
Export Report	117
Chapter 9: Loader	123
Automated Sample Loader Components	123
Loader Acquisition Controls and Settings	125
Loader Daily Workflow	127
Startup	127
Set Up Experiment	127
Acquisition	128
Loader Shutdown	132
Chapter 10: Maintenance	135
Maintenance Schedule	135
Scheduled Maintenance	135
Unscheduled Maintenance	135
SIT Flush	136
Purge Filter	136
Clean Flow Cell	137
Long Clean	137
Fluidics Shutdown	138
Cleaning the External Surfaces	139
Inspecting the Fluidics Lines	139
Replacing the Sheath Filter	139
Replacing the SIT	140
Loader Maintenance	143
Calibrating the Loader Stage and Plate	143
Cleaning the External Surfaces of the Loader	144
Loader Connections	145
Uninstalling the Loader Power Supply	145
Loader Lifting and Carrying Instructions	145

Chapter 11: Troubleshooting	147
General Troubleshooting	148
Loader Troubleshooting	150
Chapter 12: Glossary	151
Chapter 13: Specifications	155
Chapter 14: Supplies and Replacement Parts	161
Index	163

Introduction

Intended Purpose

Cytek Northern Lights™-CLC (NL-CLC) flow cytometer system (including SpectroFlo® CLC software) is intended for use with in vitro diagnostic flow cytometric applications to identify cells in suspension. It is intended to analyze human specimens in conjunction with fluorescence-labeled antibody reagents to generate information that serves as an aid to diagnosis, monitoring, or prognosis.

Cytek Automated Sample Loader (ASL) is an optional sample loading accessory for in vitro diagnostic use by adding high-throughput sample acquisition capabilities to Cytek Northern Lights™-CLC (NL-CLC) flow cytometer.

SpectroFlo® QC Beads are intended for routine performance tracking and set up of Cytek Northern Lights™-CLC (NL-CLC) flow cytometer for in vitro diagnostic use, by providing a standardized method for quality control of the instrument optics, electronics, and fluidics.

About this Guide

This guide provides information on the Cytek Northern Lights™-CLC flow cytometer, daily workflow, SpectroFlo® CLC software features, instrument maintenance, accessories, and cytometer specifications. It also includes troubleshooting tips and service information.

Safety

Safety Symbols

The following table lists symbols used throughout this guide.

Symbol	Meaning
	Caution: hazard or unsafe practice that could result in material damage, data loss, minor or severe injury, or death
	Risk of electric shock

Symbol	Meaning
	Biological risk
	Laser radiation

General Safety

- Read all safety instructions completely before using the equipment. Keep the instructions in a safe place.
- Follow all instructions when operating the instrument.
- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- Do not place any object on top of the instrument.
- Do not block any ventilation openings.
- Do not place the unit near any heat sources such as radiators, heat registers, stoves, or other devices (including amplifiers) that produce heat.
- Use only attachments/accessories specified by the manufacturer.
- Use only with the cart, stand, tripod, bracket, or table specified by the manufacturer or sold with the equipment. When a cart is used, use caution when moving the cart/equipment combination to avoid injury from tipping over.
- Unplug the instrument when it will not be used for long periods of time.
- Refer all servicing to qualified service personnel. Service is required when the unit has been damaged in any way, such as: if the power-supply cord or plug is damaged, if the unit is dropped, if liquid is spilled onto the unit or objects fall into the unit, if the unit is exposed to rain or moisture, or if the unit does not operate properly.
- Do not expose the instrument to temperatures outside the range of 15°C to 28°C.
- Concentrations of sodium hypochlorite (bleach) higher than 10%, as well as other cleaning agents, can damage the instrument. Use a 10% solution of household bleach to clean, where indicated. A 10% bleach solution is prepared by adding 1 part household bleach to 9 parts deionized water.
 - **NOTE:** Household bleach contains 5–7% sodium hypochlorite.
- Before turning on the cytometer, visually inspect all containers. Wear standard PPE laboratory attire such as protective gloves, eyewear, and lab coat.
- Purge the sheath filter if air bubbles are visible in the sheath filter, or if the plenum or sheath container have run dry.
- Fill the sheath container as needed. Never use tap water as sheath solution. Never use surfactant-based sheath solutions.
- Do not run bleach or detergent through the sheath filter. It is difficult to remove cleaning solutions from the sheath filter.
- Check the cytometer periodically for fluid leaks or crimped lines. If evidence of a leak is detected, contact Cytex Technical Support immediately. Do not attempt to repair the instrument.
- When performing Daily QC, always select the current bead lot number.
- **Notice for EU:** Users should report any serious incident related to the use of the device to the Manufacturer and National Competent Authority.
- **Notice for outside EU:** Contact your local Cytex office for any incident related to this device.

Electrical Safety

- Do not place liquids on top of the instrument. Any spill into the ventilation openings could cause electrical shock or damage to the instrument.
- Do not use this equipment near water.



WARNING: To reduce the risk of fire or electric shock, do not expose this apparatus to rain or moisture.

- Use only the power supply cord specified by the manufacturer. The power cord of the unit is equipped with a 10A three-prong power plug. Do not remove the ground pin of the power plug under any circumstances. Make sure the plug is securely plugged into the power outlet to prevent fire. If the power supply cord needs to be replaced, the cross section area of the conductor should be at least 16 AWG. This is to prevent electric fire or shock.
- Use only the fuse specified by the manufacturer. The fuse is 250 VAC, 5 A, size 5 x 20 mm.
- Protect the power cord from being tread upon or pinched, particularly at the plug and the point where it emerges from the equipment. Be sure that the power outlet is located near the equipment so that it is easily accessible.

Biological Safety

- Empty the waste container when filling the sheath container or as needed to prevent leakage. Take care to avoid damaging the fluid level sensor in the waste tank.
- Biological samples are potentially dangerous and/or life threatening. Adhere to proper handling procedures for samples and reagents. Wear standard PPE laboratory attire such as protective gloves, eyewear, and lab coat.
- Any instrument surface in contact with biological specimens can transmit potentially fatal disease. Use universal precautions when cleaning the instrument or replacing parts.

Laser Safety

The Cytex Northern Lights™-CLC cytometer is a Class 1 Laser Product and complies with the US FDA Center for Devices and Radiological Health 21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007. No laser radiation is accessible to the user during normal instrument operation.

The lasers are fully contained within the instrument and do not require any special work area safety except when service procedures are being performed. These procedures are performed only by Cytex service personnel.



CAUTION: Modification or removal of the optics covers or laser shielding could result in exposure to hazardous laser radiation. To prevent damage to skin and eyes, do not remove the optics covers or laser shielding, or attempt to service the instrument where laser warning labels are attached.

Use of controls or adjustments or performance of procedures other than those specified in this user's guide may result in hazardous radiation exposure.

To prevent exposure to laser radiation:

- Do not defeat any safety interlocks on the instrument.
- Do not use controls, make adjustments, or perform procedures other than those specified in this user's guide.
- Do not attempt to perform service procedures on the lasers.

Technical Support

For instrument support within the US, call 1-877-922-9835 (option 1).

Outside the US, call 1-510-657-0102. In Europe, call +31207653440.

Email Cytex at technicalsupport@cytekbio.com.

Visit our website, www.cytekbio.com, for up-to-date contact information.

When contacting Cytex, have the following information available:

- Serial number
- Any error messages
- Details or screen shots of recent system performance
- SpectroFlo® CLC software version and system firmware version, located in the Help module

Overview

Northern Lights™-CLC System

The Northern Lights™-CLC system consists of the Cytex Northern Lights™-CLC flow cytometer and a computer workstation running SpectroFlo® CLC software for acquisition and analysis. SpectroFlo® QC Beads are also provided for instrument QC.

The cytometer is an air-cooled, compact benchtop instrument. It is equipped with up to three lasers, 38 detection channels for fluorescence, and three channels for scatter (blue laser FSC, blue laser SSC, and violet laser SSC). Sheath and waste fluids are contained in either 4-L tanks, included with the system, or 20-L cubitainers. Software indicators notify you when the sheath is getting low or the waste is getting full. The fluidics system includes a plenum for storing sheath, allowing you to fill the sheath tank during operation.

An optional high-throughput sample loader is available to automate sample delivery and acquisition. The Automated Sample Loader (ASL) offers added compatibility for standard 96-well plates, deep 96-well plates, and 40-tube racks. See [“Loader” on page 123](#) for more information.

The workstation is a dedicated USB-compatible PC with monitor, keyboard, and mouse. It runs Microsoft® Windows® 10 or Windows® 11 Pro with a 64-bit operating system, which is required for SpectroFlo® CLC software.




Cytometer Overview

Solid-state lasers transmit light through a flow cell where particles in suspension are focused, single file for interrogation by the laser. Proprietary, high-sensitivity, semiconductor detector APD arrays are equipped with up to 16 channels per laser to capture the emission spectra of dyes that emit in the 365 to 829-nm wavelength range. The resulting fluorescence and scatter signals are then collected and converted into electronic signals. On-board electronics convert these signals into digital data that can be acquired and recorded on the workstation.

The cytometer power button is located on the left side panel in the upper right area. When the cytometer is powered on, the power button is illuminated.

The front panel opens on hinges to the left to reveal the fluidics system. The SIT door to the right of the front panel opens to reveal the sample injection tube (SIT) assembly. The top cover opens to reveal the optics plate.

Front of Cytometer




Do not place any object on top of the instrument.



Do not place liquids on top of the instrument. Fluid leaking into the cytometer could cause electrical shock or damage to the instrument.

Back of Cytometer

Allow 20.0 cm (8 in) between the back of the cytometer and the wall for proper ventilation without an air duct. Allow 10.0 cm (4 in) with an air duct.




Fluidics System

Sample Injection Port and Sample Injection Tube

Sample, contained in a standard 12 x 75-mm tube, enters the cytometer through the sample injection tube (SIT) that is contained within the sample injection port (SIP). The 12 x 75-mm sample tube snaps into place under the SIP requiring no additional tube retention support. The SIT extends from the SIP during acquisition and retracts when the cytometer is not acquiring.

When using the optional loader, a 96-well plate, 96-well deep-well plate, or 40-tube rack can be used for sample delivery in place of individual 12 x 75-mm sample tubes. See [“Loader” on page 123](#) for more information.



Fluid Containers

The cytometer draws sheath solution directly from a 20-L sheath cubitainer or the 4-L sheath tank provided by Cytek. It expels waste into an empty 20-L cubitainer or the 4-L waste tank provided by Cytek.


The included fluidics tanks are contained in a holding reservoir located on the left side of the cytometer. The 4-L tank with the transparent fluidic line is for sheath solution. The 4-L tank with the orange fluidic line is for waste.



Fluid Flow


The fluidics are driven by vacuum. An accumulator vessel is the source of vacuum for the system. Sheath solution is drawn into and stored in the sheath plenum before passing through a sheath filter, where debris and contaminants are removed. Before reaching the flow cell, the sheath stream passes through a degasser, which removes air bubbles. After passing the laser interrogation point, the combination of sheath solution and sample travels to the waste container.

Sheath and waste fluid levels are monitored by sensors. The waste level sensor is located underneath the waste tank cap. The sheath level sensor is located underneath the sheath plenum cap. Both sensors are monitored by the software.



Fluidics Components

The following figure shows the fluidics components.




The following table describes the fluidics components.

No.	Component	Description
1	Plenum pump	Pulls sheath from the sheath tank to fill the plenum
2	Vacuum pump	Maintains the vacuum in the accumulator
3	Plenum	Storage vessel for sheath fluid before it flows to the sheath filter
4	Degasser	Removes air bubbles from the sheath fluid
5	Sheath filter quick connects (x3)	Sheath filter fluid input, fluid output, and vent line quick-connects
6	Sheath filter	Filters debris and particles from the sheath fluid
7	Accumulator	Vacuum source for the fluidics system

Optics

Unlike conventional flow cytometers that direct specific bandwidths of fluorescence light into discrete detectors or photomultiplier tubes (PMTs), the Northern Lights™-CLC uses a solid-state, multi-channel, narrow-beam detector array for each laser. Each array can be configured with up to 16 detectors that are used to capture a part of the emission spectrum from each particle passing through the laser beam. The detector channels from up to three lasers are used to capture the entire emission spectra from each fluorescent-labeled particle. Spectral deconvolution (unmixing) algorithms calculate the contribution of the known individual fluorophore's spectra to the total collected signal.

For excitation, a proprietary flat-top laser design enables a constant power distribution across the width of the sample core stream.



Optical configurations are as follows:

Laser	Excitation	Channels for detection	Detector names
Violet	405 nm	16	V1–V16
Blue	488 nm	14	B1–B14
Red	640 nm	8	R1–R8

The wavelengths detected by each detector (channel) increase across the array. See the table on [page 108](#) for details.

Loader Overview

The Automated Sample Loader (ASL) is an automated loading system that mixes samples and delivers tube racks and plates to the cytometer for acquisition. The loader can be included as an

option on a new system, or it can be ordered and installed at a later time by field service. The loader offers consistent performance to meet application requirements.




The ASL resuspends samples using an orbital shaker. In addition to supporting a variety of 96-well plates, including deep-well plates, it supports a 40-tube rack. Loader settings are user-adjustable, including mixing speed, mixing duration, the number of SIT flushes between samples, and data record delay time. Loader settings are Default (standard), High throughput, and Low carryover.

Software Overview

SpectroFlo® CLC software allows you to acquire and analyze samples and adjust instrument settings. Once you log into the software, a Get started menu appears.

Get Started Menu

The Get started menu provides seven modules that allow you to perform various functions. These same seven modules are also accessible across the upper-right corner of each module screen.



The following table describes options in the Get started menu.

Module	Description
QC & Setup	Daily QC ensures that the instrument is in optimal condition for use. Run SpectroFlo® QC Beads daily to assess system performance and to adjust settings to account for day-to-day variation. Levey-Jennings reports keep track of trends in system performance. Setup allows you to create reference controls. See “QC & Setup” on page 33 for information.
Acquisition	The Acquisition module allows you to create experiments to acquire and analyze data. Experiments can be created through a guided wizard or created from previously saved templates. See “Acquisition” on page 47 for information.
Worklist module	The worklist is used for batch sample data collection and report generation of assay(s). See “Worklist” on page 111 for information.
Extra Tools	Here, FCS files can either be unmixed or compensated using virtual filters. See “Unmixing and Compensation” on page 93 for information.
Library	The library allows you to store experiment templates, worksheet templates, user settings, fluorescent tags, SpectroFlo® QC Bead information, label information, and keywords. See “Library” on page 67 for information.
Preferences	Software preferences can be changed to customize the software. Default plot sizes, fonts, gate colors, print layouts, statistics box table option, and more can all be changed in the Preferences. See “Preferences” on page 75 for information.

Module	Description
Users	The Users module contains user management options and administrative controls. See “Users” on page 86 for information.


About Experiments

The Acquisition module provides the tools necessary to acquire data, such as the acquisition controls used to start, stop, and record data, and the instrument controls used to set the threshold and adjust the detector gains. See [“Experiment Display” on page 48](#) for more information on these controls. Experiments contain the fluorescent tags and labels used in the experiment, the stopping criteria, and the groups of tubes run, which can include the reference control group. You can create groups for your samples, if you wish, to conveniently organize samples by type or staining panels, for example.

Opening an Experiment

When you click Acquisition from the Get started menu, the Acquisition Experiment menu (below) is displayed, allowing you to open a default or template experiment, create a new experiment, or import an experiment. A wizard walks you through the steps to create a new experiment.

■ **NOTE:** By default the Acquisition menu in the left pane is collapsed, showing only the icons for Experiment, Worksheet, Cytometer. To expand the menu to show the labels, click the arrows (>>) at the bottom of the pane.



Experiments can be created using several different methods. The following table describes the options in the Acquisition Experiment menu:

Method	Description
Default	Opens a new experiment with one group containing one tube and a set of labels and fluorescent tags in a default experiment worksheet template. The default experiment is user configurable. It is the quickest way to begin sample acquisition.
New	Opens the New Experiment Wizard to guide you through creating an experiment.

Method	Description
Template	Allows you to select from a list of saved experiment and assay templates (see page 24).
Import	Imports an experiment ZIP file that was exported.
My Experiments	Allows you to select from a list of saved experiments. NOTE: Original experiments can be duplicated without data, which is equivalent to opening an experiment template. Right-click an experiment in My Experiments and select Duplicate (with or without data).

Completed experiments can be accessed through the My Experiments option in the Acquisition Experiment menu. Use the column headers to sort the list of experiments. For every tube recorded, two FCS files are saved, one raw and one unmixed. Use My Experiments to open experiments you already ran, as you may want to review the data or acquire more samples. You can also export experiments from My Experiments (below). A ZIP file is exported, containing all the raw data files, if applicable (and unmixed files for unmixed experiments), as well as the worksheet templates and experiment template.

My Experiments			
 Import  Export			
Experiment	Created By	Date Created	Date Modified
Experiment_005	Admin	August 08, 2018 - 17:00 PM	August 08, 2018 - 17:00 PM
Experiment_004	Admin	August 08, 2018 - 14:47 PM	August 08, 2018 - 14:47 PM
Experiment_003	Admin	August 08, 2018 - 13:31 PM	August 08, 2018 - 13:31 PM

FCS Files


FCS files generated from an experiment are stored in the Export folder by default, or the folder you set as the default. See [“Storage Preferences” on page 82](#) for information. Experiments can contain the following types of FCS data files for each tube run:

- raw data files only (for samples that were acquired in an experiment)
- raw data files + unmixed data files (for samples that were acquired and unmixed live during acquisition)
- unmixed data files only (for samples that were unmixed post acquisition)

Experiment Templates

Use the Save As option above the experiment’s tube/group (hierarchy) list to save the current experiment as a template, which can then be used for running similar experiments. Experiment templates include fluorescent tags used in the experiment, reference controls, groups/tubes, labels, worksheets, and stopping criteria. Templates are saved in the library. To open and use a

template, select Template from the Acquisition Experiment menu. See [“Experiment Templates” on page 73](#) for more information on experiment templates.




About Worksheets

Worksheets are used to visualize the data in the experiment. Each experiment requires at least one worksheet. A worksheet allows you to view the data in plots during acquisition, as well as perform analysis functions. Worksheets contain the tools necessary to create plots, gates, annotations, statistics, and the population hierarchy. Worksheets are saved with the experiment and can be saved separately and reused across experiments.

Two types of worksheets are available—worksheets for raw data and worksheets for unmixed data. You must select the appropriate worksheet to view the corresponding type of data. When viewing raw data, the parameters on the pots in a raw worksheet reflect the channel names, for example, B1-A, R1-A, V1-A. When viewing unmixed data, the parameters on the plots in an unmixed worksheet reflect the fluorescent tags, for example, PerCP-A.

Opening a Worksheet


To select a worksheet, click **Worksheet** in the Acquisition menu. The Select a worksheet menu appears (below). You can open a new raw or new unmixed worksheet. These worksheets open with a single FSC vs SSC plot. Use the worksheet tool bar to add plots and other elements. All worksheets are saved as template files (WTML) and can be opened using the **Open Worksheet Template** option. You can also import worksheets that were exported, as well as save, print, and save a worksheet PDF.



The following table describes the options in the Acquisition Worksheet menu:

Method	Description
Default Raw	Opens a default raw worksheet that can be used for experiments where reference controls will be acquired. Do not overwrite this worksheet. Always use Save As to save this worksheet with a new name.
New Raw	Opens a new raw worksheet.
New Unmixed	Opens a new unmixed worksheet.
NxN plots unmixed	Uses the fluorescent tags selected when setting up the experiment to create a worksheet with multiple plots, displaying each fluorescent parameter against each other. This worksheet allows you to check for potential unmixing errors and correct as needed.
Open Worksheet Template	Allows you to select from a list of saved worksheet templates. A default raw and default unmixed worksheet are provided.
Import	Imports a worksheet template that was exported.
Save, Save As, Print, Save PDF	Saves the worksheet, saves the worksheet with a new name, prints the worksheet, saves a PDF of the worksheet.

You can have multiple worksheets open at a time. The currently displayed worksheet appears with a blue line under the worksheet name. Because you can select different worksheets for different groups or tubes in an experiment, each tube will have a worksheet associated with it.




Each user can define a default raw worksheet and a default unmixed worksheet. Open the default worksheet and set it up for your experiment, then select Save As to save this worksheet with a new name. The worksheet will be available to select when you create an experiment. You can use this worksheet, open a template worksheet, or create a new worksheet.

All worksheets are saved in the library. See ["Worksheet Templates" on page 72](#) for more information on worksheet templates.

About Reports


The report is a static data display sheet that can display the data from one or more tubes. Experiments can be saved as Assays. An assay needs to contain at least one report. For more information, see [“Create Reports” on page 64](#).



Startup & Shutdown

Filling the Sheath and Emptying the Waste

The color-coded sheath and waste quick-connects and the waste level sensor connector are located at the lower-left corner of the front panel.



NOTE: We recommend using the 20-L cubitainers instead of the 4-L tanks for systems with a loader.

Filling the Sheath

Fill the sheath container with manufacturer-provided sheath solution, MilliQ™ water, phosphate-buffered saline (PBS), or deionized (DI) water.

Sheath can be drawn from either the supplied 4-L sheath tank or directly from a 20-L cubitainer.

Sheath solution can be added to the sheath container while the instrument is running.



Before turning on the cytometer, visually inspect all containers for leaks or cracks. Wear the recommended protective laboratory attire such as protective gloves, eyewear, and lab coat.



Fill the sheath container as needed. Use only the appropriate sheath solution. Never use tap water or surfactant-based sheath solution.

Filling Sheath into a Cytex 4-L Sheath Tank or a 20-L Cubitainer

- 1 Remove the sheath fluidics line cap from the cubitainer or sheath tank lid from the sheath tank.
- 2 Add the appropriate sheath solution.
- 3 Replace the fluidics line cap or sheath tank lid. Do not over-tighten.
- 4 If the cytometer is powered on and the software is connected, verify that the software sheath indicator is green.

Warm Up Time Left: 00:29:31



Sheath



Waste



Cytometer



Loader

Emptying the Waste

Waste can be expelled into either the supplied 4-L waste tank or directly into an empty 20-L cubitainer.



Empty the waste container when filling the sheath container or as needed to prevent leakage. The software indicator for waste will be yellow or red when the container needs to be emptied. Take care to avoid damaging the fluid level sensor in the waste tank.



Biological samples are potentially dangerous and/or life threatening. Adhere to proper handling procedures for samples and reagents. Wear standard PPE laboratory attire such as protective gloves, eyewear, and lab coat during this procedure.



Always treat the contents of the waste container with household bleach (10% of the total volume). Contents of the waste container may contain biohazardous material.

Removing Waste from a Cytex 4-L Waste Tank or a 20-L Cubitainer

- 1 Disconnect the waste line orange quick-connect from the cubitainer cap or 4-L waste tank. Disconnect the waste level sensor.
The waste level sensor connector for the cubitainer is on the cubitainer cap. The waste level sensor connector for the 4-L tank is on the front of the cytometer.
- 2 Remove the waste cap from the cubitainer or the lid from the 4-L waste tank, taking care not to damage the liquid level sensor.
- 3 Dispose of the waste per local regulations.
- 4 Add 2 L of undiluted bleach to the waste cubitainer, or 400 mL of bleach to the waste tank.
- 5 Replace the waste cap/lid to the container. Hand-tighten the cap/lid until it is fully closed.
- 6 Reattach the waste line and level sensor line to the cap/lid and front of the cytometer.
- 7 If the cytometer is powered on and the software is connected, verify that the software waste indicator is green.

Starting Up the System

- 1 Turn on the workstation, then turn on the cytometer.


■ **NOTE:** Ensure that a tube containing 1 mL of deionized (DI) water is loaded on the SIP before launching SpectroFlo® CLC software. The tube is required for the SIT depth calibration and to flush liquid through the flow cell to remove bubbles that may have formed.

- 2 Launch SpectroFlo® CLC software and log in by entering your user name and password and clicking **SIGN IN**.

You can start typing your user name to display a list of names beginning with the letter(s) you type.

The cytometer initialization procedure begins. Sheath fluid is flushed through the fluidics lines to prevent any saline buildup, and the system calibrates the SIT depth.

- 3 Select **QC & Setup** from the Get started menu.



- 4 Check the status indicators in the lower-right corner of the screen.

- Ensure the indicator for **Connected** is a green checkmark. It may take a few minutes for the indicators to update.


Warm Up Time Left: 00:29:31 ✓ **Sheath** ✓ **Waste** ✓ **Cytometer** ✓ **Loader**

If the indicator shows the instrument is not connected, check to ensure that the USB connection between the cytometer and workstation is plugged into the appropriate ports. See [“Back of Cytometer” on page 17](#).

- Ensure the status indicators for sheath and waste are green before proceeding.

Fluid Indicator	Meaning	What to do
Yellow sheath	Sheath tank is low. Acquisition will continue for 5 minutes before the sheath is empty. You can refill the sheath during this time.	Fill the sheath tank (see “Filling the Sheath” on page 29).
Red sheath	Sheath tank is empty.	
Yellow waste	Waste tank is nearing capacity.	Empty the waste tank (see “Emptying the Waste” on page 30).
Red waste	Waste tank is full.	

- 5 Wait 30 minutes before running samples. See [“Performing Daily QC” on page 33](#). We recommend running a tube of DI water during the warm-up period. You cannot record a file until the 30 minute warm-up period is complete.



■ **NOTE:** Administrators can select Enable data collection during system warm up under Preferences > Cytometer. This bypasses restrictions imposed by default during instrument warm-up, such as data acquisition and Daily QC.

Shutting Down the System

Run the shutdown procedure at the end of each day that you use the instrument.

The shutdown procedure flushes the flow cell and sample lines with 10% bleach solution, DI water, 25–50% Contrad 70, and DI water. The procedure takes approximately 5 minutes. The SIT will remain extended from the SIP at the end of shutdown to ensure the SIT does not dry and form clogs.

To shut down the system:

- 1 Follow the Fluidics Shutdown procedure. See [“Fluidics Shutdown” on page 138](#).
- 2 When the shutdown procedure is complete, exit SpectroFlo® CLC software by clicking the X in the upper-right corner of the application window.
- 3 Turn off the cytometer and workstation.

QC & Setup

Daily QC

Run Daily QC using SpectroFlo® QC Beads prior to acquiring samples to ensure that the cytometer is performing optimally. Daily QC assesses the instrument's optical alignment and the system performance drift by measuring rCVs and gains needed to place the beads at the target locations for each detector. During QC, laser delays and area scaling factors are optimized and gain settings are adjusted to account for day-to-day instrument variability. Upon completion of Daily QC, a QC report is generated. QC reports can be reviewed under the Reports tab.

Performance can be tracked and charted over time in the Levey-Jennings tab. The software can be configured to display a warning if the QC result on the QC report exceeds user-defined criteria. See ["Alarm Ranges" on page 45](#).

Performing Daily QC

- 1 Allow 30 minutes to pass after turning on the system to ensure the optics compartment is warmed up.
- 2 Prepare SpectroFlo® QC Beads (1 drop of beads in 0.3 mL of sheath solution or PBS).

■ **NOTE:** Always prepare beads in the same solution used for the sheath solution on the instrument. The bead diluent and instrument sheath solution must match. If you use DI water to prepare the beads, the beads will begin to degrade within a few hours. Do not reuse the beads prepared in DI water.

The SpectroFlo® QC Beads are 3-µm hard-dyed, polystyrene beads that have a single fluorescence intensity. They can be excited by each laser and emit fluorescence in all detector channels.

3 Select QC & Setup from the Get started menu.



4 Select the current bead lot from the Bead Lot menu.

Each time you open a new lot number of SpectroFlo® QC Beads you must import the bead lot ID into the library so it is accessible when you run QC. Bead lot files can be downloaded from the Resources section at <http://www.cytexbio.com>.



Different bead lots have different fluorescence intensities. Always select the correct bead lot when performing Daily QC.

5 Start acquisition:

- Tube: Load a 12 x 75-mm tube of the beads onto the SIP. Select Start to begin acquisition.



- Loader: Ensure the correct plate type is selected (96 U-, V-, Flat-Bottom, deep well, or 40-tube rack). Click in the plate or tube in the tube rack to select the well or tube where the QC beads

are located. Well A1 is selected by default, but you can choose any well. Select Eject, load the plate onto the stage, then click Load followed by Start to begin acquisition.

 **Daily QC**

Daily QC: Missing

Performed on

Bead Lot: 2003

Carrier Type:  96 U-Bottom

Select QC well:

-	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Ready

Start Abort Eject

As the instrument begins acquiring the QC beads, they appear in the scatter plot. The laser delays are initially set to 0, then optimized thereafter. The performance measurements are established and compared to the pass/fail criteria (see [“Pass/Fail Criteria” on page 36](#)).

The procedure takes approximately 3 to 5 minutes to complete. Once acquisition is complete, two SIT Flushes are automatically performed to clear the beads from the sample line.

The following message is displayed when Daily QC passes. To view the QC report, click View Report (see [“QC Report” on page 36](#)).

 **Setup Completed**

Daily QC has completed and PASSED!

Close

View Report

If QC fails, remove the tube or plate and follow the guidelines that appear. The recommended solution will vary depending on the reason for the failed test.

Daily QC Failed

Daily QC failed to place the QC beads on scale.

Please try the following:

1. Make sure the concentration of your beads sample is appropriate (between 100 - 1000 events per second).
2. Go to Acquisition to run your beads sample in an experiment and make sure your sample is appropriate.
3. Remix your beads sample run QC again.

If the problem persists, please contact Cytex Service.

OK

You are now ready to run reference controls, if applicable, or acquire samples.

QC Report

At the completion of Daily QC, a QC report is generated. The report includes the following sections:

- The header section contains the Pass/Fail status of the run, name of the instrument, instrument configuration, date the Daily QC was run, user who ran the Daily QC, instrument serial number, and SpectroFlo[®] QC Beads lot and expiration date.
- The results section contains the gain, gain change, median fluorescent intensity of the QC bead, %rCV, and a pass/fail indicator for each detector channel. The center wavelength of the detector is shown in parentheses next to the detector name.
- The Laser Settings section contains the laser delays for all non-primary lasers, and area scaling factors for all lasers and the FSC detector.

Pass/Fail Criteria

The pass/fail criteria are the following:

- %rCV must not exceed 6% for the FSC channel
- %rCV must not exceed 8% for the SSC and SSC-B channels
- %rCV must not exceed 6% for the third channel of each laser (V3, B3, and R3)
- % delta gain change for all channels must not exceed 100% from the last Daily QC run performed by Cytex Service personnel.

QC reports are automatically exported as CSV files to the Setup folder (C:\CytexbioExport\Setup). The number of reports listed in the Reports screen can be set in the Preferences. See [“QC Setup Preferences” on page 84](#) for more information.

An example QC Report is shown.

Daily QC Report			
Setup Status:	PASSED	Date:	October 28, 2017 - 17:03 PM
Cytometer Name:	MyCyto	User:	Admin
Configuration:	3-Lasers-V16-B14-R8	Serial Number:	R0001
QC Beads			
Lot ID:	1002	Expiration Date:	December 31, 2019

Laser	Detector (nm)	Gain	Gain Change	Median (x1000)	% rCV	Status
Blue	FSC	174	-26	1,843.2	2.57	✓
Violet	SSC	342	9	2,087.8	4.45	✓
Violet	V1 (428)	381	55	202.2	3.96	✓
Violet	V2 (443)	212	19	205.7	3.97	✓
Violet	V3 (458)	201	17	202.5	4.11	✓
Violet	V4 (473)	153	19	244.1	3.98	✓
Violet	V5 (508)	197	13	302.0	4.04	✓
Violet	V6 (528)	248	12	243.1	4.04	✓
Violet	V7 (549)	233	13	182.7	4.01	✓
Violet	V8 (571)	256	19	123.1	3.79	✓
Violet	V9 (594)	251	14	102.5	3.80	✓
Violet	V10 (618)	381	18	91.0	3.78	✓
Violet	V11 (664)	638	46	72.6	3.73	✓
Violet	V12 (692)	974	59	60.6	3.71	✓
Violet	V13 (720)	530	36	31.5	3.75	✓
Violet	V14 (750)	531	40	21.0	3.88	✓
Violet	V15 (780)	793	75	10.8	5.96	✓
Violet	V16 (812)	461	35	4.2	7.93	✓
Blue	B1 (508)	231	2	14.1	2.62	✓
Blue	B2 (528)	242	-10	38.7	2.03	✓
Blue	B3 (549)	221	-10	94.4	1.63	✓
Blue	B4 (571)	240	-13	134.2	1.66	✓

Red	R7 (783)	895	-9	79.4	5.85	✓
Red	R8 (812)	326	0	39.3	6.11	✓

Laser Settings

Laser	Laser Delay	Area Scaling Factor
Violet	-24.95	1.19
Blue	0.00	1.20
Red	27.50	0.85

FSC Area Scaling Factor: 1.24

Specifications

FSC	% rCV:	< 6	(Recommended)
SSC	% rCV:	< 8	(Recommended)
V3	% rCV:	< 6	(Recommended)
B3	% rCV:	< 6	(Recommended)
R3	% rCV:	< 6	(Recommended)
All Channels	% Gain Change:	< 100	(Recommended)

Reference Controls

Reference controls, obtained by running single-stained and unstained samples, provide the individual fluorescence spectra necessary to unmix the data. Either beads or cells can be stained for use as reference controls. These controls can be acquired in the reference group of the experiment during acquisition, or they can be acquired as reference controls in the QC & Setup module. If reference controls are acquired in the QC & Setup module, they are stored and can be used as reference controls for unmixing in subsequent experiments.

A wizard guides you through recording reference controls.


Running Reference Controls

To create reference controls you will need to select the fluorescent tags, choose the control type (beads or cells), then label the fluorescent tags (for example, using CD nomenclature). You can also append existing FCS files when running reference controls.

Run Daily QC to ensure the instrument is performing optimally before running reference controls.

- 1 Select Reference Controls from the QC & Setup module.
- 2 Select New Reference Controls from the Reference Controls tab.

A wizard opens allowing you to create new reference controls. You can also append existing FCS files.



- 3 Select the fluorescent tags. The left pane displays the fluorescent tag groups found in the library.
 - Click the arrow to the left of the fluorescent tag group name to view the fluorescent tags in the group. (The default fluorescent tag groups are Blue Laser, Red Laser, Violet Laser, Fluorescent Proteins, and Viability Dyes and contain a list of commonly used fluorescent tags excited by their respective laser).
 - From the expanded list of fluorescent tags, select the fluorescent tags used in the experiment. Once selected, the fluorescent tags appear in the selection pane on the right. You can select fluorescent tags by dragging and dropping, double-clicking, or using the Add button. Multiple tags can be chosen at one time. Confirm the tags selected, then click Next.

■ **NOTE:** The list of fluorescent tags can be edited in the library. You can use the library to add fluorescent tags that are not present in the default list. See “Fluorescent Tags” on page 67 for more information.

4 Define the unstained control(s).

Choose either single-stained fluorescent tags, or new, separate, unstained for the control tubes. If using a control that has a positive population only, you will need an unstained control of the same type as the stained.

Either beads or cells can be stained and defined as control types. If you select Use new, separate, unstained tube(s), you can add more unstained tubes, as necessary. Then select the specific unstained tube in the Negative Control column of the Fluorescent Tags table below.

- 5 Define the control type (beads or cells) for each fluorescent tag in the Fluorescent Tags table.
- 6 (Optional) Enter labels associated with the fluorescent tag for identification and tracking.
- 7 If applicable, enter the lot number(s) of the reference controls.

■ **NOTE:** If you selected Label/Lot Specific Unmixing in the Acquisition Preferences (see [page 76](#)), the software will search the library and experiment reference groups for reference controls that have the same fluorescent tag, label, and lot information in order to use the corresponding control for unmixing.

Define Unstained Control(s)

☐ Use from Single Stained Fluorescent Tag Tubes
☒ Use new, separate, unstained tube(s)

Name	Control Type
Unstained_1	Beads
Unstained_2	Beads

+ Add X Remove

Fluorescent Tags

Fluorescent Tag	Control Type	Label	Lot
FITC	Beads		
PE	Beads		
PerCP	Beads		

8 Once the controls have been defined, click Next.

9 Load the control sample onto the SIP and click Start to preview the sample data. This step allows you to ensure that all populations for all control samples are on scale. If necessary, adjust gain settings.

Gain setting for all channels for each laser detector module can be adjusted. Adjust the detector using the all channels % feature. Adjust the same percentage for all lasers to maintain the fluorochrome signature. To save this customized setting, click the Save button in the Instrument control pane and rename.

Instrument Control

User Settings: CytekAssaySetting (Cytek) Save Save As

GAIN THRESHOLD SIGNAL LASERS

FSC SSC SSC-B

50 50 50

Violet Blue Red

V1 V2 V3 V4 V5

4 11 13 13 12

V6 V7 V8 V9 V10

9 9 5 4 3

V11 V12 V13 V14 V15

1 0 0 0 0

V16

0

All Channels %: 0

■ **NOTE:** CytekAssaySetting was generated by measuring the optimal resolution of human lymphocytes stained with CD8 and CD4 labeled with various fluorochromes. These settings provide a place to start for most immunophenotyping applications.

The Adjust Settings screen allows you to view the data to ensure that the positively stained fluorescent particles are on scale. Adjust the threshold and FSC and SSC gains. FSC gain can be adjusted from 1-1,000. SSC and detector channel gains can be adjusted from 10-10,000. If the positive population is off-scale for any detector channels, lower the gain setting for that detector module, and then lower by the same percentage for each laser in the instrument. This will preserve the fluorochrome signature.

If the positive population is not sufficiently separated from the negative population within a specific channel, verify that CytekAssaySetting is selected, as these settings ensure adequate resolution. If the positive population is still not clearly separated from the negative population, review the sample preparation and staining procedure.

Once gain settings have been confirmed, unstained and reference controls are ready for acquisition.

■ **NOTE:** Dim markers may not separate from the negative population regardless of how much the gain is increased.



10 Preview the data for the remaining controls by clicking Stop, then:

- Tube mode: Loading the next control onto the SIP and clicking Start.
- Plate mode: Selecting the next control well and clicking Start.

11 Select Next when you are satisfied with the gain settings. Proceed to running controls.

12 If you are running in Tube mode, place a tube of the appropriate single-stained particles on the SIP.

13 Click Record to begin acquiring.

Make sure to follow the order listed in the left-hand panel.

During acquisition the spectra plot for each fluorescent control is displayed. The plots show all the channels across all lasers in the x-axis vs mean fluorescence intensity (MFI) of the fluorescent tag.

If the acquisition criteria is not met within 15 minutes, the run will stop.



14 During acquisition, obtain spectral information by moving the polygon gate on the FSC-A vs SSC-A plot to include the population of interest.

Hold down the Ctrl key while adjusting the gate to move the polygon gates for all the scatter plots at once. The gated population appears in the histogram, which is set approximately to the peak emission channel of the fluorescent tag to be acquired. The emission spectrum of the population is displayed in the spectrum plot.

Adjust the positive gate on the histogram. The software automatically displays the emission spectrum of the positive particle in the spectrum plot. SpectroFlo® CLC software sets the default gate near or on the peak emission channel. The gate can be selected manually. It is best to set the gate on the brightest emission, as this can make distinguishing the positive and negative populations easier and provides better visualization of the spectrum.

■ **NOTE:** Unmixing results are unaffected by the position of the interval gate in the spectrum plot.

Readjust the positive and/or negative gate on the histogram, if necessary.

15 Continue recording each control.

16 Select Save to save the reference controls.


FCS files for saved reference controls are stored in the export Setup folder by default, or the folder you set as the default. See ["Storage Preferences" on page 82](#) for information.

Setting Reference Controls as Benchmarks for Reference Control QC

You can select reference controls to be used as benchmarks for QC. The reference controls you select as benchmarks allow you to check the quality of your reference controls prior to unmixing. If a reference control is suboptimal, the benchmark will indicate whether the reference control spectra appears as expected. If you choose not to designate reference controls as benchmarks, the reference controls run during QC & Setup or during the experiment will be used as-is without a quality control check.


NOTE: Before saving a reference control as a benchmark, ensure that it is fully optimized.

- 1 Select Reference Controls from the QC & Setup module.
- 2 Select Benchmark from the Reference Controls tab.



- 3 Select the reference controls that will be the benchmarks.

Click All to select all the controls in the list.



- 4 Click Save.

Updating Reference Controls


You may wish to update the reference controls if any of the following occur:

- Major service performed on the instrument
- Fluorochrome exhibiting signs of instability resulting in changes in its emission spectrum
- Staining protocols change, for example different buffers or fixatives are used

The Reference Controls tab displays the reference controls saved in the library. Click the arrow next to the control name to display the details.

To Update Reference Controls

- 1 Select Update Reference Controls from the Reference Controls tab in the QC & Setup module.
A wizard opens allowing you to update reference controls.



- 2 Follow steps 3 through 11 in “Running Reference Controls” on page 38.


Levey-Jennings Tracking

Levey-Jennings tracks the gain and CV's for all detector channels and the laser delays and area scaling factors for each laser over time. This allows you to view the system's performance. Select the parameter(s) you wish to track.

The graphs in the report show you random errors or shifts and trends in the data for each parameter. Data from the last 30 days, 3 months, or 12 months can be included in the reports.

- 1 Select a parameter checkbox to display the Levey-Jennings plot.
- 2 Use Save Selection to save the LJ tracking settings.

3 To export Levey-Jennings data to a .csv file, click Export Data.



Gain Settings

The amount of signal amplification applied to each detector channel can be modulated by increasing or decreasing the amount of gain applied. The gains for every detector channel can be saved and are collectively known as the user settings. User gain settings are stored as a ratio against the Daily QC. Every time Daily QC is performed, User Settings will be adjusted accordingly.

Alarm Ranges

You can set an alarm to warn you when the gain and %rCV exceeds the passing criteria that you define. This changes the outliers (shown in red) in the LJ graphs. Select Alarm Range from the Cytometer QC tab, then adjust the SD range (plus or minus) for individual detectors for each laser.

QC & Setup Cytometer QC				
QC & Setup Acquisition Extra Tools Library Preferences Users Help				
Cytometer QC Alarm Range Save				
Parameter	Alarm Criteria	Min	Max	
Gain	+/- 3SD	0	0	
Violet	+/- 3SD	0	0	
Blue	+/- 3SD	0	0	
Red	+/- 3SD	0	0	
% rCV	+/- 3SD	0.00	0.00	
Violet	+/- 3SD	0.00	0.00	
Blue	+/- 3SD	0.00	0.00	
Red	+/- 3SD	0.00	0.00	
Laser Delay	+/- 3SD	0.00	0.00	
Violet	+/- 3SD	0.00	0.00	
Red	+/- 3SD	0.00	0.00	
Area Scaling Factor	+/- 3SD	0.10	3.00	

Acquisition

Raw vs Unmixed Data

SpectroFlo® CLC software saves flow cytometry data in the FCS 3.1 format. Data is saved in both raw and unmixed formats. Raw data contains all the fluorescence information from each detector (ie, V1, V2, V3, etc). Each detector channel is designated by its excitation laser and position in the array. For example, B3 is the third channel of the blue laser detector array.

Unmixed data contains all the fluorescence information from each fluorescent tag in the experiment. To unmix data, single stained controls (or reference controls) for each of the fluorescent tags (as well as an unstained control) are required. During unmixing a mathematical algorithm is used for the decomposition of the fluorescent components in the sample using the reference controls. Parameters in unmixed data will display as the fluorescent tag name along with their associated labels.

The Acquisition module provides the tools that allow you to create an experiment. An experiment is a set of tubes, instrument settings, acquisition criteria (stopping rule), fluorescent tags, labels, and worksheets designed for the acquisition of samples. See [“About Experiments” on page 23](#).

New and saved experiments can be created or accessed in the Experiments tab of the Acquisition module.

Unmixing and Compensation

Raw FCS files can be spectrally unmixed in the following ways:

- Reference group from the experiment – Reference controls collected as FCS files within the experiment can be used to unmix using the Unmixing wizard in the Acquisition module.
- Reference controls run in QC & Setup module – Reference controls run in the QC & Setup module can be used to unmix using the Unmixing wizard in the Acquisition module.
- Unmixing from the Extra Tools module – FCS files collected from different experiments can be unmixed in the Extra Tools module. FCS files can be imported and unmixed in this module.

Raw FCS files can also be compensated with the conventional method using the Virtual Filters tab in the Extra Tools module. Detector channels can be binned together to simulate the analysis of the data as if it were acquired using a filter. See [“Virtual Filters” on page 106](#) for more information.


Setting Up an Experiment

An experiment can be saved as a template or created for one time use. Setting up the experiment in SpectroFlo® CLC software involves:

- 1 (Optional) Providing a name and description for the experiment. A default name is provided.
- 2 Specifying the fluorescent tags used in the experiment.
- 3 Defining the reference group with associated reference tags, labels, and lot numbers, as needed.
- 4 Adding labels and lot numbers to each of the fluorescent tags.
- 5 Adding custom keywords. Custom keywords can be defined in the Library.
- 6 Selecting an acquisition worksheet—either new or a template.
- 7 Defining acquisition criteria (stopping rule based on events, time, or volume).

Acquisition Experiment Overview

The Acquisition module provides the necessary elements for data collection within the experiment. Click the Experiment icon in the far left pane to open a template, the default, or a new experiment using the wizard.




Experiment Display

The experiment display in the Acquisition module includes the following panes. To show, hide, or undock (float) these panes from the experiment panel, click the corresponding icons in the top-right corner of the pane.

Group-Tube List and Hierarchy


The samples are listed in the upper left of the screen. Samples can be organized into groups. Use the (+) Tube and (+) Group icons to add tubes and groups. If running from a loader, click Add Plate or Rack to add a new sample carrier. Click Plate View to display a graphic image of the carrier instead of the list of groups.

Click **Save** to save changes to the experiment, or click **Save As** to save an experiment template. Click **Edit** to edit the experiment. Worksheets can be applied to the experiment, groups, or individual tubes.



Acquisition Control

The Acquisition Control pane allows you to start, stop, and pause acquisition, record data, and restart acquisition counters. The acquisition controls are enabled when a tube is present on the SIP. To show, hide, or undock (float) this pane from the experiment panel, use the dock/undock and hide icons in the top-right corner.



The following table describes the controls in the Acquisition Control pane.

No.	Control	Description
1	Start/Record/Pause/Stop/Restart	Start and Record are enabled when a tube is present on the SIP. Select Start to start acquisition. Select Record to record data. Record can also start acquisition. Select Pause to pause recording. While paused you can adjust the flow rate. Select Record again to continue. Select Stop to stop acquisition. Select Restart to restart the acquisition counters. All events and results displayed are refreshed. Stop and Restart are enabled once Start is selected. Stop and Pause are enabled once Record is selected.
2	SIT Flush	Select to perform a SIT Flush

No.	Control	Description
3	Flow Rate	Select Low (15 µL/min), Medium (30 µL/min), or High (60 µL/min). The exact flow rate is displayed.
4	Event Rate, Abort Rate, Threshold Count, Time Elapsed	Displays the real-time counts during acquisition.
5	Events to Display	Enter the number of events to display during acquisition.

Instrument Control

The Instrument Control pane consists of the Gain, Threshold, Signal, and Lasers tabs for use in adjusting the instrument settings.

User Settings allows you to select CytekAssaySetting, Default, or any saved user settings for the experiment. We recommend using CytekAssaySetting as a starting point. This setting provides the optimal resolution for each channel, accommodates bright signals, and minimizes spread. While using CytekAssaySetting, you will need to only adjust FSC, SSC gains, and Threshold.

Instrument Control

User Settings: **CytekAssaySetting**

GAIN THRESHOLD SIGNAL LASERS

Save Save As

FSC SSC SSC-B

299 656 456

Violet Blue Red

V1 V2 V3 V4

1,153 406 335 243

V5 V6 V7 V8

225 254 283 224

V9 V10 V11 V12

259 561 548 735

V13 V14 V15 V16

491 544 844 322

All Channels %: 0

THRESHOLD SIGNAL LASERS

Threshold Operator: ☐ Or ☒ And

Channel Threshold

FSC 100,000

GAIN THRESHOLD SIGNAL **LASERS**

Window Extension: 2.00

FSC Area Scaling Factor: 0.87

Laser Area Scaling Factor

Violet 1.19

Blue 0.93

Red 0.99

Laser Laser Delay

Violet -24.33

Blue 0.00

Red 28.75

SIGNAL LASERS

Scatter Channels

FSC ☒ Area ☒ Height ☐ Width

SSC ☒ Area ☒ Height ☐ Width

SSC-B ☒ Area ☒ Height ☐ Width

Fluorescence Channels


Laser	Area	Height	Width
Violet	☒	☒	☐ V9
Blue	☒	☒	☐ B8
Red	☒	☒	☐ R5

The following table describes the tabs in the Instrument Control pane.

No.	Tab	Description
1	Gain	Gains can be adjusted for all detector channels for all lasers using the gain spinboxes. FSC gain can be adjusted from 1–1,000. SSC and fluorescence detector gains can be adjusted from 10–10,000. To change the value that the gain increments, see “Acquisition Preferences” on page 75 . Use All Channels % to increase/decrease all gains for a selected laser by the percentage you select.
2	Threshold	Use the Threshold tab to set the threshold parameter and minimum threshold channel value. Multiple parameters can be set as a threshold using either the AND or OR operator. Use OR when at least one parameter is available.
3	Signal	Use the Signal tab to select area, height, or width for each signal. Area and height can be selected for all channels. Width can be selected for only one channel per laser.
4	Lasers	Use the Lasers tab to set the area scaling factor and laser delay. These values are automatically set and updated in all user settings upon completion of the Daily QC. ■ NOTE: If you run large cells and the lowest FSC gain setting is not low enough to see your cells, lower the FSC area scaling factor (for example, 0.5).

Worksheet

The worksheet allows you to view the data in plots and create plots, statistics, population hierarchy, and gates. Click the worksheet icon in the far left pane to select, import, save, and print worksheets.



Worksheet Toolbar

A toolbar at the top of the worksheet area allows you to undo/redo, zoom; create plots, gates, statistics, population hierarchy, annotations; and save, print, and save a PDF of the worksheet. Hover the cursor over an icon to see a description and keyboard shortcut.



Plots

Four plot types can be created in the worksheet:

- dot plots
- pseudocolor plots (density plots)
- histogram plots
- spectral plots

To change the properties of a plot, right-click the plot and select Properties. You can select the plot type, parameters, scale, background color, and labels.

Plot Properties

General

Plot Gate: Tube_001_All Events ▾

Plot Type: Pseudocolor Plot ▾

Parameters

X Axis Parameter: FSC-A ▾

X Axis Scale: Linear ▾

Y Axis Parameter: SSC-A ▾

Y Axis Scale: Linear ▾

Layout

Width: 300 Height: 300

Density Plot Options

Density Levels: 15 ▴ ▾

Miscellaneous

Background Color:

☐ Include Plate Name

☐ Include Group Name

☒ Include Tube Name


☒ Include Population Name

☐ Include Custom Title

Gates

Gate types include:

- Rectangle
- Ellipse
- Polygon
- Interval
- Hinged Quadrants
- Staged Quadrants (select and drag an offset handle to move the quadrant segment up or down)
- Binary



Click and drag handle to move quadrant segment up or down.

Gate Properties

Gate properties can be changed by right-clicking the gate. You can change the name of the gate, the color, and gate boundary line weight. You can also select whether to display the count and/or the % parent events within the gate, as well as the gate parameters.

Gate Properties

Gate Name:

P1

Gate Color:

☐ Count

☐ % Parent

Gate Boundary Line Weight:

Normal

Gate Borderline Color:

Parameters

X Axis Parameter

FSC-A

Y Axis Parameter

SSC-A

Boolean (Logical) Gates

Select multiple gates and right-click to open a menu to:

- Intersect the gates with the AND operator – events that are present in all of the selected gates are part of the intersected gate population.
- Join the gates with the OR operator – events that are present in at least one of the gates are part of the joined gate population. You can assign a new color for the joined gate or keep the original gate color (see “Gate Properties” on page 54).



Statistics

To create a statistics box, click the Statistics icon in the worksheet toolbar, then click in the worksheet area.



Select the population checkbox next to the populations that have stats to display. To add a statistic, select the statistic from the Statistics Variable list.

Select the parameter you would like to add for the statistics. Multiple parameters can be selected at once.

The software offers a counts/ μ L statistic that can be calculated for any gate.

To adjust the precision of the statistics, select the decimal place in the Decimal Places table. To remove a statistic, right-click the column header and select Delete.

You can add keywords to statistics.

Create Statistics Table

Statistics

Header

Select a Control to View

Available Controls

Custom keywords

Standard keywords

Keyword	☒ Select All
\$PAR	☐
\$TOT	☐
\$MODE	☐
\$NEXTDATA	☐
\$BYTEORD	☐
\$DATATYPE	☐
\$BEGINDATA	☐

Once the statistics are complete, you can export the stats as a CSV file to the location you choose. Statistics output can be appended or overwritten if you export the stats into the same CSV file.

You can add keywords to the statistic table for exporting. Click the Header tab in the Create Statistics Table window and select the keywords to export.

Create Statistics Table

Statistics

Header

Title

Statistics Table 1

Add Statistics

☒ Population

☒ Select All

☒ Count

☒ % Parent

☒ % Grand Parent

☒ All Events

☒ Single cells

☒ Lymphocytes

☒ CD3+

Statistics Variable

Parameter

Decimal Places

Median

rSD

% rCV

Mean

Max

Min

FSC-H

SSC-H

SSC-B-H

FITC-A

PE-A

PerCP-Cy5.5-A

Median

rSD

% rCV

Mean

Max

Min

0

0

2

0

0

0

Add Statistics


Population	Count	% Parent	% Grand Parent	% Total	Count/uL
All Events					
Single cells					
Lymphocytes					
CD3+					
CD4-CD8+					
CD4+CD8+					
Q1-R3					
CD4+CD8-					
CD3+CD8+					
CD3+CD4+					
CD3-					
CD(16+56)+					
Q2-R2					

Cancel

Save

Next


To export statistics from multiple tubes, select a group or multiple tubes, then right-click to select Batch Analysis. PDF and CSV files with data for multiple tubes can be exported.




Creating a Default Experiment

Use the default experiment option to quickly begin an experiment. The default experiment contains a set of fluorescent tags, however, all parameters are user configurable.

- 1 Click Default in the Acquisition Experiment menu.




The default experiment opens with a default raw worksheet, plate/40 tube rack.



- 2 Choose carrier type.
- 3 Click Edit.

The same wizard appears as when creating a new experiment. Follow the steps to select the fluorescent tags, add groups, and define makers, keywords, and acquisition settings/worksheets.


See “Creating a New Experiment” on page 58 for details on defining the information in the experiment wizard.




Creating a New Experiment

Selecting New in the Experiment menu opens the New Experiment wizard. The wizard walks you through the steps to create a new experiment.

- 1 Click New in the Acquisition Experiment menu.



- 2 The Create New Experiment wizard opens. Specify a name for the experiment or use the default name. (Optional) Type in a description.



- 3 Click the arrow to the left of the group name (laser) in the Library pane on the left to display the list of its fluorescent tags. Select the fluorescent tags used in the experiment and click (+) Add to

add them to the Selection list on the right. You can also double-click the tag to add it to the Selection list.

To quickly find a fluorescent tag, type the tag name in the Type to filter text box. A default list of fluorescent tags for each group is available in the library. See “Fluorescent Tags” on page 67.


You must select all fluorescent tags present in the experiment, as this will determine which reference controls are to be used during spectral unmixing.

The screenshot shows the 'Create New Experiment' window with the 'Fluorescent Tags' tab selected. The 'Name' field contains 'Experiment_005'. Below it is a 'Description (optional)' field. A 'Library' list on the left contains various fluorescent tags: PE, PE-CF594, PE-Dazzle594, PE-eFluor 610, PE-Texas Red, PE-Alexa Fluor 610, PE-Cy5, PE-Cy5.5, PerCP, PerCP-Cy5.5, BB700, PE-Alexa Fluor 700, and PerCP-eFluor 710. A 'Type to filter' search box is positioned above the library list. To the right of the library list are 'Add' and 'Remove' buttons. On the right side of the window, a 'Selection' list shows 'FITC', 'PE', and 'PerCP', with 'PerCP' highlighted. Above the 'Selection' list is a '3 Tags Added' status and a 'Clear All' button. At the bottom left are 'Cancel' and 'X' buttons, and at the bottom right is a 'Next' button with a right arrow.

To remove individual fluorescent tags from the Selection list, click to select the tag, then click Remove. To remove all fluorescent tags, click Clear All.

- 4 Once all fluorescent tags have been chosen from the Library list, confirm the list in the selection pane, then click Next.
- 5 Ensure the correct Carrier Type (Manual Tube), plate, 40 tube rack, or plate type is selected, then create groups for your samples by selecting  Group. Add tubes to the groups, plate, or 40 tube rack.

- 6 Create groups for your samples by selecting  Group. Add tubes to the groups.




- 7 Select  Reference Group if you are intending to unmix with all or some controls acquired in this experiment.

This creates a list of reference control tubes for each fluorescent tag specified as part of the experiment.

■ **NOTE:** If you plan to unmix the samples using only reference controls run in QC & Setup, step 7 is not necessary.

■ **NOTE:** When using controls from the reference library and running controls live, you will delete the controls to be used from the library and run the other controls. When unmixing, you will add the controls from the library to unmix all the controls.

- 8 **IMPORTANT:** Define an unstained control for autofluorescence by selecting its control type (beads or cells). The unstained control needs to be of the same type and prepared in the same way as the samples, as this will ensure accurate unmixing and autofluorescence quantitation.



Fluorescent Tags				
Fluorescent Tag	Control Type	Label	Unstained	
FITC	Cells			
PE	Cells			
PerCP	Cells			

- 9 If applicable, select Define Extra Unstained Control(s) for Spillover Calculation to use a different unstained control to calculate spillover for your reference controls. Then enter a name and control type for this extra unstained control.

For example, if test samples are cells and the reference controls are beads, all with only positive peaks, you will need to run a separate tube of negative beads for the spillover calculation. An extra unstained control is not needed if your unstained autofluorescence control (and sample) is the same type as the reference controls.

Create Reference Group

Define Unstained Control(s) for Autofluorescence Extraction

Name: Unstained Control Type: Cells

☒ Define Extra Unstained Control(s) for Spillover Calculation

Name	Control Type
Unstained CompBeads	Beads

- 10 Select the control type (beads or cells) for the single-stained reference controls.
- 11 (Optional) Enter the label (for example, CD nomenclature) that is conjugated to the fluorescent tag.
- 12 If applicable, enter the lot number(s) of the reference controls.

■ **NOTE:** If you selected Label/Lot Specific Unmixing in the Acquisition Preferences (see [page 76](#)), the software will search the library and experiment reference groups for reference controls that have the same fluorescent tag, label, and lot information in order to use the corresponding control for unmixing.

Create Reference Group

Define Unstained Control(s) for Autofluorescence Extraction

Name: Unstained Control Type: Cells


☐ Define Extra Unstained Control(s) for Spillover Calculation

Fluorescent Tags					
Fluorescent Tag	Control Type	Label	Lot	Unstained	
FITC	Cells				
PE	Cells				
PerCP	Cells				

■ **NOTE:** Use the red trash can icon to delete an individual tube from the reference group. This may be necessary if you wish to mix and match references acquired in this experiment with reference controls run in QC & Setup. Any stored controls you plan to use should be deleted from the reference group.

13 Click Save.


Once the reference group has been created, entries for each of the references will be displayed. Each of the reference group tubes will have an icon (with the letter R) associated with it under the reference group.



14 If necessary, continue adding tubes, click Next when all tubes are created.

15 Add markers/labels to the remaining sample tubes before continuing. They can be chosen from the Labels list on the right, typed directly into the table, or copied and pasted. Labels can be added at the group or tube level and can be applied to multiple cells selected at once. Labels are required for reference controls if you selected Label/Lot Specific Unmixing in the Acquisition Preferences (see [page 76](#)). Click Next when all the tubes are labeled.

■ **NOTE:** If you selected Label/Lot Specific Unmixing, the software will search the library and experiment reference groups for reference controls that have the same fluorescent tag, label, and lot information in order to use the corresponding control for unmixing.



16 (Optional) Enter custom keywords and click Next.

Custom keywords can be added at the experiment, group, or tube level. You must define the custom keywords in the Library before you can add them to an experiment (see ["Keywords" on](#)

page 69 for information). Drag and drop the keywords from the Keywords list on the right to the experiment, group, or tube and enter keywords values as needed. You can also copy and paste custom keywords across different tubes in this wizard.

Edit Experiment

Fluorescent Tags Groups Markers **Keywords** Acquisition

Value

Name	Keyword	Value	Keyword	Value	Keyword	Value	Keyword
▼ Test-Keywords							
▶ Reference Group							
▶ Setup							
▼ Samples							
WB1	Sample ID	A01	Patient Name	ABC	Age	25	
WB2	Sample ID	A02	Patient Name	XYZ	Age	46	

Keywords

- ▶ Subject ID
- ▶ LymphSet 6C
 - Sample ID
 - Patient Name
 - Age

17 Select the acquisition settings and worksheet(s). The worksheet menu lists all the worksheets for the given user.

- Select the Default Raw Worksheet (Raw) for the Reference Group and for the sample groups if you plan to perform post-acquisition unmixing.
- Select the Default Unmixed Worksheet (Unmixed) or any user-created unmixed worksheet for your sample groups if you are performing live unmixing. Worksheets can be selected at the experiment, group, or tube level.

Select the stopping gate, storage gate, number of events to record, Stopping volume (μL), Stopping Criteria, Stopping Time (seconds). Acquisition stops when the first of the stopping criteria is met (time, volume, number of events). These criteria can be selected at the experiment, group, or individual level.

- If acquiring beads, we recommend collecting 5,000 singlet events.
- If acquiring cells, we recommend collecting 10,000 to 20,000 events of the desired population.

■ **NOTE:** The number of events to acquire depends on the target population. For example, you may need to acquire 10,000 to 20,000 events to get 2,000 of the desired population. Approximately 1,000 to 2,000 events is needed in both the negative and positive populations of each control for accurate unmixing.

- 18 Select Tube/Well Specific User Settings if you want to apply specific user settings to individual tubes/well. If left unchecked, all tubes/wells in the experiment will use the same experiment user settings.

Edit Experiment

Fluorescent Tags Groups Markers Keywords Acquisition

☒ Tube/Well Specific User Setting

	Worksheet	Stopping Gate	Storage Gate	Events To Record	Stopping Volume (ul)	Stopping Criteria	Stopping Time (sec)	User Setting
int_001	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000	3,000	☐ Count & Volume	10,000	CytekAssaySetting (Cytek) ▼
up_001	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000	3,000	☐ Count & Volume	10,000	CytekAssaySetting (Cytek) ▼
ie_001	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000	3,000	☐ Count & Volume	10,000	CytekAssaySetting (Cytek) ▼

Previous Cancel Save and Open

- 19 Once the worksheet and stopping criteria have been defined, click **Save and Open** to open the new experiment.

To make any changes to the experiment, click **Edit** above the group/tube hierarchy.

- 20 To acquire controls and samples and perform live unmixing, see [“Live Unmixing”](#) on page 96.

Creating an Assay

Create New Assays

Assay components

An Assay is similar to an experiment template, but is formatted to be used in a Worklist for batch sample data collection. An assay includes the following elements: tube(s) (one or more tubes), fluorescence tags, marker labels, worksheet(s) (optional), report(s) (at least one), preview time, data collection stopping rules, Carrier type, and Keywords

Create Reports

Experiments can be saved as Assays. An assay needs to contain at least one report.

- 1 Click the plus sign to the right of **REPORTS** on the experiment toolbar to create a report.


2 Double-click the corresponding Report to edit the report name.



3 Click on report.

You can edit the report and add elements such as plots, gates, statistics, header and footer information. Adding worksheet elements to the report is done by copying from the worksheet desired and pasting into the report.

4 Copy the desired plots and gating hierarchy from the worksheet and paste in the report.



NOTE: Report vs. Worksheet

- The Report is a static data display sheet that can display the data from one or more tubes.
- The Worksheet is a dynamic data display sheet, which displays the data in the current tube with tube point pointed to it and can only display the data from one tube at a time.

5 Add Header, Footer, and Signature in Report.

- Select the header, footer, and signature icon buttons on the toolbar to add headers, footers and signature boxes to the report.
- Enter relevant information in the header and footer as needed, such as the name of the report, company logo, etc.




- 6 Right-click in the header or footer on the Report to select the information from a dropdown list to display, such as Sample ID, Operator, Data recording date, etc.

7 Color TBNK Assay				
Institution: Cylekbio Cytometer SN: 16898 Software: SpectroFlo® CLC & 1.0.1 (07142021)		Assay: 7C B Operator: Admin Page: Page 1 Of 1		Sample ID: Approve date: Status:
Gate	% Parent	% Grand Parent	Count/uL	Absolute Count
7C B_Lymphocytes	15.71	13.66	59.44	619.22
7C B_CD3+	62.59	9.83	37.21	386.95
7C B_CD3+CD4+	84.41	52.84	31.41	326.65
7C B_CD3+CD8+	13.77	8.62	5.12	53.29
7C B_CD19+	45.82	17.14	10.19	105.97
7C B_CD(16+56)+	52.58	19.67	11.69	121.60
7C B_CD20+	45.61	17.06	10.14	105.48
7C B_CD19+ AND CD20+	45.61	17.06	10.14	105.48
7C B_CD19+ OR CD20+	45.82	17.14	10.19	105.97

Save an experiment as an assay

Users with proper privileges can use Save As-Assay to convert an experiment into an assay template (Assay), which can be used in the Worklist.



Library, Preferences, and Users

Library

The library contains information for various elements used for the experiments. Information saved in the library includes SpectroFlo[®] QC Beads lots, fluorescent tags, labels, user settings, worksheet templates, experiment templates, and keywords. Information stored in the library can be saved, exported, and imported for reuse.

QC Beads

SpectroFlo[®] QC Bead lot IDs and expiration dates can be imported, exported, or removed from the library. The QC bead lot for the beads used for Daily QC must be saved in the library. Select QC Beads to see a list of QC bead lots.

Library QC Beads							
QC Beads Fluorescent Tags Labels User Settings Worksheet Templates Experiment Templates Keywords	QC Beads Import Export Remove <table> <tr> <th>Lot ID</th><th>Expiration Date</th></tr> <tr> <td>1002</td><td>December 31, 2021</td></tr> <tr> <td>1001</td><td>December 31, 2020</td></tr> </table>	Lot ID	Expiration Date	1002	December 31, 2021	1001	December 31, 2020
Lot ID	Expiration Date						
1002	December 31, 2021						
1001	December 31, 2020						

Fluorescent Tags

Fluorescent tags are the designation given to each distinct fluorescent molecule that can be detected by the system. This includes for example, fluorophores, fluorescent proteins, and fluorescent viability dyes. Each unique fluorophore run on the instrument must be given a fluorescent tag name.

By default, several groups of fluorescent tags are pre-installed with the software—Blue Laser, Red Laser, Violet Laser, Fluorescent Proteins, and Viability Dyes. These groups contain the most

commonly used fluorophores excited by the system's on-board lasers. Additional tags can be added to these groups.

The screenshot shows the 'Library' window in the SpectroFlo CLC software. The left sidebar contains a navigation menu with options: QC Beads, Fluorescent Tags (selected), Labels, Keywords, Spillovers, User Settings, Loader Settings, Worksheet Templates, Experiment Templates, and Backup & Restore. The main area is titled 'Fluorescent Tag Groups' and includes buttons for Import, Export, New, Edit, and Delete. Below this is a table of existing groups. To the right, the 'Current Group: Blue Laser' is selected, showing a table of fluorescent tags with columns for Name, Emission Wavelength, Laser Excitation, and Display Name.

Name	Created By	Date Created	Date Modified	Name	Emission Wavelength	Laser Excitation	Display Name
Blue Laser	System	August 11, 2021 - 14:29 PM	August 11, 2021	BB660	667	488	BB660
Red Laser	System	August 11, 2021 - 14:29 PM	August 11, 2021	PerCP	677	488	PerCP
Violet Laser	System	August 11, 2021 - 14:29 PM	August 11, 2021	cFluor BVG680	680	488	cFluor BVG680
Fluorescent Proteins	System	August 11, 2021 - 14:29 PM	August 11, 2021	BB700	693	488	BB700
Viability	System	August 11, 2021 - 14:29 PM	August 11, 2021	PE-Cy5.5	695	488	PE-Cy5.5
				cFluor B690	690	488	cFluor B690
				PerCP-Cy5.5	695	488	PerCP-Cy5.5
				PE-Fire 700	695	488	PE-Fire 700
				PerCP-Vio700	704	488	PerCP-Vio700
				cFluor BVG710	710	488	cFluor BVG710

You can create groups or individual fluorescent tags by selecting Add. Both groups and individual fluorescent tags can be imported or exported.

■ **NOTE:** Selecting New will create a new group while selecting Add will allow you to add an individual fluorescent tag to a group.

To edit the properties of a fluorescent tag that you added, select the fluorescent tag of interest and select Edit. Properties that can be edited include fluorescent tag name, laser excitation wavelength, emission wavelength, and display name. The default tags that are included with the software cannot be edited or deleted. However, they can be exported and imported across different systems running SpectroFlo® CLC software.

If the fluorophore is known by another name or identified by a different spelling, those additional names or spellings can be added as synonyms.

Any groups that you created can also be edited in this window. The default groups cannot be edited or deleted.

Labels

Fluorescent tags can be conjugated or attached to proteins that can specifically bind to other proteins on the cell surface or within the cytoplasm. They can also be inherently fluorescent, such as fluorescent proteins that can be fused to a variety of cellular proteins using molecular cloning techniques. The proteins that are either bound or attached to fluorescent tags can be designated as labels. The software comes with an initial set of pre-installed labels that are categorized as CD

Markers, Chemokines, Chemokine Receptors, and Cytokines. Additional labels can be added by using  Add in the right pane.

Library

Labels

QC Beads

Fluorescent Tags

Labels

User Settings

Worksheet Templates

Experiment Templates

Keywords

Label Groups

Import

Export

New

Edit

Remove

Add

Name

Created By

Date Created

Date Modified

Label

CD Markers

System

July 13, 2018 - 10:43 AM

July 13, 2018 - 10:43 AM

CD1a

Chemokine

System

July 13, 2018 - 10:43 AM

July 13, 2018 - 10:43 AM

CD1b

Chemokine Receptors

System

July 13, 2018 - 10:43 AM

July 13, 2018 - 10:43 AM

CD1c

Cytokines

System


July 13, 2018 - 10:43 AM

July 13, 2018 - 10:43 AM

CD1d

CD1e

New label groups and individual labels can be created by clicking New. Label groups can also be imported and exported for use on other systems. The default labels cannot be edited or deleted.



 New Label Group

Group Name:

Description (optional):

Library

Type to filter

- CD Markers
 - CD158J
 - CD158K
 - CD175
 - CD175s
 - CD176

 Add

 Remove

Selection

Drop labels here!

 Clear All

 Cancel  Save


Keywords

Keywords allows you to type, import, or export user-defined, custom FCS keywords by group. These keywords can be assigned in an experiment at the experiment, group, or tube level and are exported with the FCS files.

The following types of keywords can be created:

- Numeric
- String
- Boolean
- Selectable numeric

- Selectable string



Standard Keywords

The following table describes the standard FCS 3.1 keywords found in SpectroFlo® CLC FCS files.

Keyword	Description
\$FIL	Tube name
\$VOL	Total sample volume (microliters)
\$DATE	Date that the FCS data set was acquired
\$BTIM	Beginning time of data acquisition
\$ETIM	End time of data acquisition
\$OP	Operator name
\$INST	Institution name
\$CYT	Name of cytometer (Northern Lights)
\$CYTSN	Cytometer serial number
\$SPILLOVER	Spillover matrix

Custom Keywords

The following table describes the custom keywords found in SpectroFlo® CLC FCS files. Custom keywords are not required as part of the FCS 3.1 standard.

Keyword	Description
APPLY COMPENSATION	True or False to indicate whether compensation is enabled when FCS file is loaded into SpectroFlo® CLC.
CHARSET	Encoding. UTF8.
CREATOR	Software name and version
FSC ASF	FSC area scaling factor
GROUPNAME	Name of sample group
LASER{x}NAME	Name for each laser
LASER{x}DELAY	Laser delay for each laser

Keyword	Description
LASER{x}ASF	Area scaling factor for each laser
P{x}DISPLAY	Preferred display scale
THRESHOLD	Threshold channel and value
TUBENAME	Name of sample tube
WINDOW EXTENSION	Window extension used to record the tube
USERSETTINGNAME	User setting name used to record the tube

Users can change the value of a custom keyword, then select Auto Update for that keyword to update the keyword in subsequent experiments.

The screenshot shows the 'Library Keywords' interface. The left sidebar has a 'Keywords' tab selected. The main panel displays a table of keyword groups. The 'Default Keywords' group is expanded, showing a list of keywords. The 'Auto Update' checkbox for the 'Age' keyword is highlighted with a red box.

Name	Suffix	Type	Value	Auto Update	Description
Age		Numeric	0	☒	
Sample Type		String		☐	
Gender		Selectable String	M	☐	
Dilution Factor		Numeric	1.00	☐	

Spillovers

Spillovers shows the spillover matrices saved during experiments. These saved matrices can be used for other experiments.

The screenshot shows the 'Library Spillovers' interface. The left sidebar has a 'Spillovers' tab selected. The main panel displays a table with columns: Name, Created By, Date Created, and Description. The 'Spillover Matrix' section is visible on the right.

Name	Created By	Date Created	Description

User Settings

User Settings are the set of gain settings, threshold, and signal type for all detector channels. The date when it was created, as well as the name of user who created it are also saved. The settings are saved from the Instrument Control Pane in the Acquisition module. The name and description can be modified in this tab. Administrators can check the Shared checkbox to allow all user to see

and use the settings and worksheets created by other users. These settings can be saved from the Acquisition screen.

If you wish to reset the system-defined Default user settings to their original settings, select the Default user setting, then click Restore Default.

Library User Settings

QC Beads, Fluorescent Tags, Labels, **User Settings**, Worksheet Templates

User Settings

Import Export Remove Restore Default

Name	Created By	Date Created	Description
CytekAssaySetting	Admin	August 15, 2018 - 10:14 AM	Cytek Assay Setting
Default	Admin	August 15, 2018 - 10:14 AM	Default

Loader Settings

Loader settings can be saved and are stored in the library. The settings name, date the settings were created, version, as well as the name of user who created them are displayed. Click Description to add information about the settings.

Loader settings include mix time and speed, sample return, number of SIT flushes, and record data delay time.

Library Loader Settings

QC Beads, Fluorescent Tags, Labels, Keywords, Spillovers, User Settings, **Loader Settings**, Worksheet Templates

Import Export Delete

Name	Created By	Date Created	Version	Description
------	------------	--------------	---------	-------------

Worksheet Templates

All worksheets created in the Acquisition module are saved in the library and can be accessed through the Worksheet Templates tab. Worksheets can be exported as WTML files and imported for re-use. For more information on worksheets, see ["About Worksheets" on page 26](#).

Select a user(s) to view the worksheet templates created by that user(s). Click the pen tool to edit the description of the selected worksheet. Administrators can check the Shared checkbox to allow all users to see and use the worksheet templates created by other users. SuperUsers can view, import, and export worksheet templates, but cannot delete templates created by other users.

To view a worksheet, select it, then click View. To remove a worksheet that is no longer needed, select it, then click Remove.

You can also set the default raw and unmixed worksheets back to their original defaults by selecting the worksheet, then clicking Restore Default.

Experiment Templates

Experiment templates contain the key elements of an experiment without the data. They can be saved and stored in the library for future use. Saved Assays are displayed in the assay tab. The name, creation date and time, description, and creator information is displayed. Experiment templates can also be imported and exported from this tab. Assays can be used to create experiments. For more information on experiments, see [“About Experiments” on page 23](#).

Select a user(s) to view the experiment templates created by that user(s). Click the pen tool to edit the description of the selected worksheet. Administrators can check the Shared checkbox to allow all users to see and use the experiment templates created by other users. SuperUsers can view, import, and export experiment templates, but cannot delete templates created by other users. To delete an experiment template that is no longer needed, select it, then click Delete. The default experiment template cannot be deleted.

Backup & Restore

Users and SuperUsers can back up selected data from their accounts directly from SpectroFlo® CLC software. This includes experiments, experiment templates, worksheets, user settings, loader settings, daily QC results, and reference controls. Administrators have the additional ability to back up selected data from the different User and SuperUser accounts.

- 1 To back up data, click Backup in the right pane.

- 2 Choose a timeframe (day, week, month, or custom timeframe) that includes the data, select the user(s).

The screenshot displays the 'Library Backup & Restore' interface. The sidebar on the left contains the following items: QC Beads, Fluorescent Tags, Labels, Keywords, Spillovers, User Settings, Loader Settings, Worksheet Templates, Experiment Templates, and **Backup & Restore** (highlighted with a red box). The main area features a 'Backup' button, a 'Restore' button, and a table of users. The table has columns for 'User', 'Full Name', and 'Role'. The 'User' column includes checkboxes for selection. The right sidebar contains checkboxes for backup options: Data, Experiments, Experiment Templates/Assays, Worklists, Worksheets, User Settings, Loader Settings, Daily QC Results, and Reference Controls. The 'Day' timeframe is selected, and the start and end dates are both 3/2/2021.

User	Full Name	Role
☒ Admin	Administrator	Administrator
☐ Operator	Op	User
☐ Janelle S	Janelle S	User
☐ Alex Z	Alex Z	User
☐ Patty T	Patty T	User
☐ JXi	January Xi	SuperUser
☐ CytekAdmin	Cytek Administrator	Administrator

Preferences

The Preferences module allows you to change various functionality and display elements of the software user interface. The following section describes the options that can be changed in the Preferences module. Each section within the Preferences can be restored to its default settings by selecting Restore Default Preferences.



Restore Default Preferences

Acquisition Preferences

In the Acquisition tab, you can change the number of events displayed on plots during acquisition, the number of SIT flushes, the preview time before recording begins, the value that the gain settings increment when you click the spin-box arrow during acquisition, and the SIT Lift Distance. You can also enter default prefixes for tube, group, and experiment names. A label/lot specific unmixing feature allows you to choose specific reference controls for unmixing. The Batch Analysis options allow you to choose the output file type and tube/well results contained in the PDF.

Preferences

Acquisition

Acquisition

Worksheet

Plot

Gates

Statistics

Fonts

Annotation

Notifications

Storage

QC & Setup

Cytometer

Users

Global Information

Number of Events to Display on Plots:

2,000

Number of Automatic SIT Flush After Tube Unload:

1

Data Review Time(Sec) After Tube Recorded:

0

Record Data Delay Time(Sec) for Manual Tube:

0

Gain Spinbox Up/Down Increment (Ctrl key Held)

5

Experiment

Default Tube Name Start With:

Tube_

Default Group Name Start With:

Group_

Default Plate Name Start With:

Plate_

Default Tube Rack Name Start With:

TubeRack_

Default Experiment Name Start With:

Experiment_

Default User Setting:

CytekAssaySetting (Cytek)

☒ Ref Controls Ordered by Control Type

Unmixing Feature

☐ Label/Lot Specific Unmixing

Preferred SSC Channel for Unmixing:

SSC

Batch Analysis

Default Output File Type:

☒ PDF☒ CSV

Default PDF File Option:


☒ All in one file☐ Each tube/well has its own file

The following table describes the options in the Acquisition preferences.

Item	Description
Number of Events to Display on Plots	The number of events displayed in the pseudocolor plots, dot plots, and histograms. The default is 2,000 events.
Number of Automatic SIT Flush After Tube Unload	Choose the number (0–2) of SIT Flushes that will be performed once a tube is removed from the SIP.
Data Review Time (Sec) After Tube Recorded	The number of seconds that elapse before the tube pointer moves to the next tube after the current tube is finished recording.
Record Data Delay Time (Sec) for Manual Tube	The number of seconds before it begins recording after you click Record.
Gain Spinbox Up/Down Increment (Ctrl key Held)	Increments the gain for each detector channel by the amount indicated when you hold the Ctrl key and select the up and down arrows of the Gain spinbox.
Experiment	Enter the default name for tubes, groups, plates, tube racks (ASL only), and experiments. Edit the default name and/or append to the existing names. Select the default user settings for experiments. Ref Controls Ordered by Control Type allows you to order reference controls by type (beads or cells). If unchecked, reference controls are ordered by excitation laser, then by peak emission.
Unmixing Feature	Label/Lot Specific Unmixing allows you to use lot-specific reference controls for unmixing when appropriate. If this option is selected, the software will search the library and experiment reference groups for reference controls that have the same fluorochrome, label, and lot information in order to use the corresponding control for unmixing. Select the preferred SSC channel for unmixing (SSC for the signal from the violet laser (default), or SSC-B for the signal from the blue laser).
Batch Analysis	Select the default output file type for batch analysis reports (PDF, CSV, or both). Choose to create a single PDF for all tubes/wells, or individual PDF files for each tube/well.

Worksheet Preferences

The Worksheet tab allows you to change the way elements are displayed in the worksheet. Header and footer properties are also adjusted in the Worksheet tab.




The following table describes the options in the Worksheet preferences.

Item	Description
Population Hierarchy Window Size	Sets the default height and width for the population hierarchy experiment element.
Statistics Table Window Size	Sets the default height and width for the statistics table.
Grid	Display Grid – Toggles on/off the display of grid lines in the worksheet. Display Page Line – Toggles on/off the page break line in the worksheet. Grid Size – Modifies the size of the grid squares. Options include 1", 1/2", 1/4", and 1/8". Snap to Grid – Toggles on/off the ability for the worksheet elements to snap to and line up with the grid lines on the worksheet.

Item	Description
Page Setup	Title – Select a worksheet title from a list of default titles. The title is shown when the worksheet is printed or exported as a PDF. Show & Print Page Number – Toggles whether the page number is shown and printed. Print Header & Footer – Toggles whether the header and footer are printed. Allows you to select the text that is displayed in the left and right headers and footers. Print Grid – Toggles whether the grid is printed. Can also be set to use Page Setting. Page Layout – Toggles between landscape and portrait. Page Size – Sets the page size according to standard paper sizes. Page Margin – Sets the margins of the page to Narrow, Normal, or Wide.

Plot Preferences

The display properties of pseudocolor, dot, and histograms plots can be adjusted in the Plot tab.




The following table describes the options in the Plot preferences.

Item	Description
Default 2D Plot Size	Set the default height and width of the pseudocolor plots and dot plots in pixels.
Default Histogram Size	Set the default height and width of histograms in pixels.
Default Spectrum Size	Set the default height and width of spectrum plots in pixels.

Item	Description
Default Background Color	Set the default background color for all plots. Click to select from the available colors.
Default Plot Title	Customize the title of all plots to include the plate name, group name, tube name, population name, and/or a custom name.
Pseudocolor Plot Density Levels	Increase or decrease the number of density levels displayed in the pseudocolor plot.
Histogram Smooth	Set whether histogram distributions are smoothed.
Histogram Filled	Set whether histogram distributions are filled.
Histogram Y Axis	Set the scale of the histogram y-axis to count or percentage.
Image Export Format	Select the format (PNG, TIFF, or JPEG) and resolution (96, 200, 300, or 500 dpi) for exported images.

Gates Preferences

Gate properties can be adjusted in the Gates tab.



The following table describes the options in the Gates preferences.


Item	Description
Default Name Location	Select where the gate name is displayed with respect to the gate itself.
Show % of Parent together with Gate Name	Toggles on/off the display of the % of Parent with the gate name.

Item	Description
Show Count together with Gate Name	Toggles on/off the display of the population count with the gate name.
Gate Boundary Line Weight	Sets the thickness of the line drawn by the gate.
Default Colors for First X Gates	Set the number of gates (1–10) that will follow the color scheme detailed in the gate color table. The order in which the colors appear can be changed.
Interval Gate Default Color	Toggles on/off whether the population captured by the interval gate has a default color.
Quadrant Gate Default Color	Select whether the population captured by the quadrant gate has a default color, and if so if each quadrant has a different color.
Binary Gate Default Color	Select whether the population captured by the binary gate has a default color.

Statistics Preferences


The default degree of precision (number of decimal places) of the statistics displayed in the worksheet can be modified in the Statistics tab.

The precision for the following statistics can be adjusted: Mean, rSD, % rCV, Mean, Max, Min, SD, % CV, % Total, % Parent, and % Grand Parent.



Fonts Preferences

The Fonts tab allows you to change the font properties of each display element.




The following table describes the options in the Fonts preferences.

Item	Description
Font Locations	Select which display element's font to modify.
Text Settings	Font Family – Select the font family. Font Size – Select the font size. Color – Select the font color. Font Style – Toggles between normal and italic. Font Weight – Select normal, bold, or semibold.
Sample Text	Shows a preview of the text with the properties set in the Text Settings.

Notifications Preferences

The Notifications tab allows you to change certain notification settings in the Acquisition and Extra Tools modules.




The following table describes the options in the Notifications preferences.

Item	Description
Acquisition Module	Toggles whether to display the Save Changes pop-up window when closing an Experiment, Worksheet, and User Settings, as well as whether to display a confirmation message when canceling the experiment wizard.
Extra Tools Module	Toggles whether to display the instructional dialog boxes in the Extra Tools module.
Window Sizes and Locations	Select Restore Default to set pop-up window sizes and locations to the defaults.

Storage Preferences

The Storage tab allows administrators to set the default storage locations for the experiment FCS files and setup FCS files. This option is available only for administrators.

■ **NOTE:** Do not move files stored in these folders.




The following table describes the options in the Storage preferences.


Item	Description
FCS Files Folder	Select the folder where FCS files are saved. FCS files from experiments where only raw data was acquired, as well as FCS files from experiments where live unmixing was performed are stored here. For experiments where live unmixing is performed, both raw and unmixed FCS files are saved. - You can select to label FCS files with group names as the prefix. Note that you can set the default group name in Acquisition preferences. - You can select to organize FCS file folders by date. If not selected, FCS files are ordered by experiment name, without a parent folder indicating the date. If you change the default folder where FCS files are saved and want to move existing experiments to the specified folder, click Migrate Now.
Setup FCS Files Folder	Select the folder where FCS files generated by QC & Setup procedures are saved.

QC Setup Preferences

The QC Setup tab allows you to select the days/months of QC reports to display in the Reports section of the Cytometer QC menu in the QC & Setup module.



Cytometer Preferences



The following table describes the options in the Cytometer preferences.

Item	Description
Default Carrier type	Allows you to set a carrier as default
Fluidics	Use checkboxes to set preferences for fluidics boost and bubble/clog detection.
Cytometer Warmup	Click checkbox to enable data collection upon warmup.
Sample Carrier Configuration	You can have default, high throughput, or low carryover selected for each carrier type


Users

User accounts can be managed in the Users module. User account information and use time are stored in the module.

There are two types of user accounts—SuperUser and User. Only administrators can manage user accounts.

Managing Users

Administrators can add, remove, edit, and disable user accounts from the User tab. User passwords can also be changed/reset. The User tab lists all users and displays the type and status of each.




The screenshot shows the 'Users' module interface. On the left is a sidebar with navigation options: 'User' (selected), 'Use Time', 'User Role', and 'User Policy'. The main area displays the 'User' tab with a table of users. Above the table are buttons for 'Add New', 'Edit', and 'Delete'. The table has columns for 'User', 'Full Name', 'Role', 'Status', 'Email', 'Phone Number', and 'Password'. There are three users listed: 'Admin' (Administrator, Enabled), 'User' (User, Enabled), and 'Super User' (Super User, Enabled). Each user has a 'Change Password' or 'Reset Password' button next to the password field.

User	Full Name	Role	Status	Email	Phone Number	Password
Admin	Administrator	Administrator	Enabled			Change Password
User	User	User	Enabled			Reset Password
Super User	Super User	User	Enabled			Reset Password

Adding a New User Account

- 1 Click  Add New in the User tab of the Users module.

This option is available only for administrators.



The 'Create User' form is shown. It has a blue header with a pencil icon and the text 'Create User'. Below the header are several input fields: 'User Account ID' (with the value 'Cytek'), 'Password' (with masked characters), 'Confirm Password' (with masked characters), and 'Full Name'. Below these fields are two sections: 'User Type' with radio buttons for 'Administrator' and 'Operator' (selected), and 'Account Status' with radio buttons for 'Enabled' (selected) and 'Disabled'.


- 2 Enter a user account ID and password, then enter the password again to confirm.

The user account ID appears in the user name filed in the User's tab.

- 3 (Optional) Enter the user's full name.
- 4 Select the user type—Administrator, SuperUser, or User.
- 5 Select the user status—enabled or disabled.
- 6 Click Save.

Editing a User Account

- 1 Select the user from the User tab of the Users module, then click Edit.




- 2 You can edit or add a user name. You can also change the user type and/or account status.
- 3 Click Save.

Resetting a User Password

An administrator can change their or another administrator's password or reset an operator's password.

- 1 Select the user from the User tab of the Users module.
- 2 Click Reset Password.



User	Full Name	Role	Status	Email	Phone Number	Password
Admin	Administrator	Administrator	Enabled			Change Password
User	User	User	Enabled			Reset Password
Super User	Supr User	User	Enabled			Reset Password

- 3 A dialog opens displaying the new password.
- 4 Make note of the password, and click OK.

Removing a User Account

To remove a user account, select the user from the User tab of the Users module, then click Remove. When you remove a user account all the account information and that data acquired by that user is deleted. The FCS files saved by the user are retained.

Use Time

Select the Use Time tab of the Users module.

Use Time shows the total daily, weekly, or monthly use time (duration) and number of sessions that each user is on the system. Choose Admin, User, or SuperUser to view and export the use times for the respective user account.

User	Full Name	Day	Duration	Session Count
Super User	Supr User	2021-08-12	20s	1
Admin	Administrator	2021-08-12	36s	2
User	User	2021-08-12	25s	1
Admin	Administrator	2021-08-11	24h 59m 24s	2

To see the individual login sessions for each day, select Login Sessions. Choose Admin, User, or SuperUser to view and export the login sessions for the respective user account. The log on and log off times for each session, as well as the session duration are displayed. Administrators can manage the list of login sessions by deleting all sessions before a specified date. To delete sessions, click Manage Use Time, select the date, and click Delete.

User	Full Name	Day	Duration	Session Count
Super User	Supr User	2021-08-12	20s	1
Admin	Administrator	2021-08-12	36s	2
User	User	2021-08-12	25s	1
Admin	Administrator	2021-08-11	24h 59m 24s	2

Administrators can export daily, weekly, and monthly use times and login session to .csv files by selecting the user(s) and session dates and clicking Export or selecting Admin, User, or SuperUser and clicking Export All to export all use times or login sessions.

User Role

Administrators can edit the privileges for User and SuperUser user roles. Administrators can also create and edit new custom user roles. Administrator privileges cannot be changed.

Editing a User Role

- 1 Select User Role from the Users tab, then select the user name (role) you wish to edit.
The default user roles (and names) are Administrator, SuperUser, and User. The Administrator role cannot be edited.
- 2 Click Edit to edit the selected user role.
- 3 Edit the user role privileges, then click Save.

See “User Role Privileges” in the following section.

- 4 If you wish to restore user role privileges to the defaults, click Restore Default in the Edit User Role window.

User Role Privileges

The following table shows the privileges for the three default user types (Administrator, SuperUser, and User). Administrators can edit the privileges for the SuperUser and User. Administrator privileges cannot be edited.

Privileges	User	SuperUser	Administrator
Instrument Control			
Change fluorescence gains		✓	✓
Customize parameters (signal types)		✓	✓
Window extension		✓	✓
Area scaling factors	✓	✓	✓
Change laser delay		✓	✓
Enable data acquisition during warm up (cytometer preference)		✓	✓
Library Documents			
Manage QC bead lot		✓	✓
Create fluorescent tags	✓	✓	✓
Create label	✓	✓	✓
Create reference controls	✓	✓	✓

Privileges	User	SuperUser	Administrator
User Management			
Manage users		✓	✓
Manage all user types			✓
Create user role			✓

Creating a User Role

Administrators can create new user roles based on the privileges of the three default roles.

- 1 Select User Role from the Users tab.
- 2 Click New.
- 3 Enter a name for the new user role.
- 4 Select the base type for new role. Base types are Administrator, SuperUser, and User.

The new role can have the same privileges as the selected base type. You can remove some privileges but you cannot grant more privileges than the base type has. For example, if you want the new user role to manage QC bead lot information, select SuperUser as the base type, as the User base type does not have this privilege.

- 5 (Optional) Enter a description.
- 6 Click Save.

The new user role appears in the Name list.

Deleting a User Role


Administrators can delete a custom user role.

- 1 Select User Role from the Users tab.
- 2 Click Delete.
- 3 Click Delete to confirm deletion.

User Policy

Administrators can configure general password and login security policies.

- 1 Select User Policy from the Users tab to display the user policies.



Users User Policy

User

Use Time

User Role

User Policy

Password

Expires After: **Never** day(s)

Can Be Reused After: **0** day(s)

Password Strength: ☐ Strong ☐ Medium ☒ Weak

Login

Maximum Login Attempts: **Never**

The user account will automatically be locked after the set number of attempts. An Admin will be required to unlock the account.

If Inactive, Automatically Logout After: **Never** minute(s)

☐ Change Default Password When Login Next Time

☒ Enable Login Username Typing Prompt

- 2 Select the following:
 - if and when passwords expire
 - if and when passwords can be reused
 - password strength
 - strong is 8–20 characters with upper-case and lower-case letters, at least one number, and at least one special character
 - medium is 6–20 characters with letters and numbers
 - weak has no restrictions
 - if and when a maximum number of login attempts fail (users are locked out of SpectroFlo after the selected number of unsuccessful login attempts is reached)
 - if and when users are automatically logged out if inactive for the selected time (1 minute to 45 hours, in various increments)
 - change default password when login next time (if the default password Rainbow is set, the user will be prompted to create a new password)
 - enable login username typing prompt (as you begin to type when logging in, a list of user names with matching characters appears)

Unmixing and Compensation

Spectral Unmixing

Spectral unmixing is an important concept to understand how data is generated and analyzed using the Cytex Northern Lights™-CLC flow cytometer with SpectroFlo® CLC software. Spectral unmixing is used to identify the fluorescence signal for each fluorophore used in a given experiment.


Understanding Full Spectrum Flow Cytometry

Because fluorophores emit light over a range of wavelengths, optical filters are typically used to limit the range of frequencies measured by a given detector. However, when two or more fluorophores are used, the overlap in wavelength ranges often makes it impossible for optical filters to isolate light from a given fluorophore. As a result, light emitted from one fluorophore appears in a non-primary detector (a detector intended for another fluorophore). This is referred to as spillover. In conventional flow cytometry, spillover can be corrected by using a mathematical calculation called compensation. Single-stained controls must be acquired to calculate the amount of spillover into each of the non-primary detectors.

The ability to measure a fluorophore's full emission spectra allows the system to use a different method for isolating the desired signal from the unwanted signal. The key to differentiate the various fluorophores is for those to have distinct patterns or signatures across the full spectrum. Because the system is looking at the full range of emission of a given fluorophore, and not only the peak emission, two dyes with similar emission but different spectral signatures can be distinguished from each other. The mathematical method to differentiate the signals from multiple fluorophores/dyes is called spectral unmixing and results in an unmixing matrix that is applied to the data. While not mathematically identical to conventional compensation, the overall principal is the same. Just as for compensation, single-stained controls, identified in SpectroFlo® CLC software as reference controls, are still necessary, as they provide the full fluorescence spectra information needed to perform spectral unmixing.

One advantage of unmixing is the ability to extract the autofluorescence of a sample and treat it as a separate parameter. This is especially useful when running assays with particles that have high autofluorescence and for which that high background has an impact in the resolution of the fluorescent signals. Per experiment, you can define one unstained control or multiple unstained controls (one per group), depending on whether the multicolor samples have the same or different autofluorescence signatures.

In the following figure on the left, spectrum plots from conventional spectrum viewer shows heavy overlap between Qdot 705 and BV711 peak emission spectra. On the right, spectrum plots from Cytel Northern Lights™-CLC show distinct signatures for Qdot 705 and BV711.



Unmixing Workflows

Unmixing Overview

There are three unmixing workflows available in SpectroFlo® CLC software—two in the Acquisition module and one in the Extra Tools module:

- live unmixing during acquisition
- post-acquisition unmixing (in the Acquisition module)
- post-acquisition unmixing (in the Extra Tools module)


When data is acquired with live unmixing, references are acquired as raw data either in the experiment as part of the reference group or previously acquired in the QC & Setup module as reference controls. References for *all* fluorescent tags used in a given experiment must be present in the system in order for live unmixing of multicolor samples to occur. The live unmixing functionality allows you to visualize unmixed data during acquisition.

Reference Controls for Unmixing

Depending on when you unmix the data, you will use the following controls for unmixing.

Unmixing	Reference Controls
Live unmixing	Reference controls run in the experiment or reference controls run in QC & Setup
Post-acquisition unmixing in Acquisition module	Reference controls run in the experiment or reference controls run in QC & Setup
Post-acquisition unmixing in Extra Tools module	Any FCS files from samples run in any experiment

Multicolor samples can be acquired as raw data and unmixed post acquisition as well. This can be done in either the Acquisition module or the Extra Tools module.




Negative/Unstained Controls

In addition to positive reference controls needed for spectral unmixing, an unstained control is also necessary to assess **autofluorescence**. The unstained control needs to be of the same type and prepared in the same way as the samples, as this will ensure accurate unmixing and autofluorescence extraction, if desired. Ideally, your reference controls, negative control, and samples will all be the same sample type and prepared in the same way.

In addition to assessing autofluorescence, **fluorescence spillover** must also be determined. To correct for spillover, the unstained autofluorescence control can be used if it matches the sample and reference control type. However, if your reference controls do not match your sample type and

do not contain a negative population in each tube (have only positive peaks), you must use a separate spillover unstained control that matches your reference control type.

Controls



Live Unmixing


Samples can be unmixed during acquisition. Live unmixing can be performed with the reference group acquired during the experiment, the reference controls (run during QC & Setup and stored in the system), or a combination of both.

For each sample tube that is live unmixed, two FCS files are generated, one that is composed of raw data and one that is composed of unmixed data.

Live unmixed data can be analyzed in unmixed worksheets in the Acquisition module. Unmixed worksheets are different from raw worksheets, as they only display fluorescence information categorized into the defined fluorescent tags for each of the experiments.

To Perform Live Unmixing

- 1 Create a new experiment with fluorescent tags defined. Create a reference group in the experiment with the fluorescent tags, if there are any that have not already been stored as reference controls. See [“Creating a New Experiment”](#) on page 58 for details.



- 2 To view the data for the reference control tubes, make sure CytekAssaySetting is selected, then click Start. If necessary, use the Instrument Controls to adjust the settings so that all events are on scale. View all the controls, as well as the multi-color tube, and make any instrument adjustments to ensure populations are on scale before you begin recording.

To edit the acquisition criteria, click **Edit** at the experiment level and select the **Acquisition** tab. Or, to edit the properties of a single tube, right-click a tube and select **Tube Properties**.

■ **NOTE:** Keep in mind the more events you acquire, the longer it takes to unmix the data.

- 3 Click **Record** when you are ready to begin acquisition. Acquisition stops when the first stopping criterion is met.

■ **NOTE:** If necessary, you can pause to change the flow rate.

- 4 When all reference controls are acquired, click **Unmix** in the upper-left toolbar.
- 5 For the unstained controls, we recommend selecting **Use Control from Experiment** if unmixing with controls you acquired in the experiment.
- 6 If necessary, for the stained controls, select **Use Control from Library** if unmixing with reference controls run in QC & Setup.

Checkmarks appear for those controls coming from QC & Setup. The checkbox is only active if reference controls for those fluorescent tags are already saved with the reference controls from the QC & Setup module.

Unmix Experiment

Select Controls Identify Positive/Negative Populations QC Controls

UNSTAINED CONTROLS

☐ Use Control from Library

☒ Use Control from Experiment

Reference Group - Unstained (Cells)

Name	Control Type
Reference Group - Unstained (Cells)	Cells


STAINED CONTROLS

☐ From Library	Fluorescent Tag	Control	Unstained	Generic
☒	FITC	FITC (Cells)		☒

- 7 Click **Next**.
- 8 Use the **Identify Positive/Negative Populations** tab to include the positive and negative populations for each fluorescent tag in the appropriate gate.


Only the data plots for the samples you acquired are displayed, not for reference controls that you chose to use from the library.

■ **NOTE:** If you need to set the FSC and/or SSC axis to a log scale, select the Log checkbox.



- Move the polygon gate in the FSC vs SSC plot on the left to include the singlet population. Hold down Ctrl to move all the polygon gates at once.
- Move the positive interval gate in the histogram to include the positively stained population. Move the negative interval gate to include the negative population.
- Move the interval gate on the spectrum plot on the right to select the channel that exhibits the brightest fluorescence intensity. This channel is the peak emission channel for the fluorescent tag.

■ **NOTE:** If one of the controls is questionable or does not contain sufficient data, you can reacquire it or append to it, then unmix again.



- (Optional) To see how the reference controls run in the experiment compare to the benchmark reference controls, click Next.

■ **NOTE:** For information on creating benchmarks, see [“Setting Reference Controls as Benchmarks for Reference Control QC”](#) on page 43.


Two options allow you to view how the two reference controls compare—Spectral Profile and Similarity™ Matrix.

- Spectral Profile displays the emission spectrum of the unmixing controls against benchmark spectra designated by the user. The benchmark reference control spectra appear in red and the reference controls appear in black.

A similarity index appears to the right of the plots. A value closer to 1 indicates similar spectral signatures, while a value closer to 0 indicates spectral signatures that are not similar. For comparison to benchmarks, the value should be close to 1. If the value is below 0.97, it will be flagged with a yellow warning symbol. This indicates a mismatch of the unmixing control spectra with the benchmark spectra. If the similarity index falls below this value, it is imperative to check the unmixing control against the reference spectra provided in the fluorochrome guideline found in the Help menu. See “[Similarity™ Matrix](#)” in the following section.

The similarity index can also be used when viewing and comparing the full spectral signatures of any two dyes. Dyes with similar spectral signatures can be challenging to resolve. In this case, it is best to use dyes with a similarity index ≤ 0.98 . See spectrum.cytekbio.com for a full spectrum viewer tool.

If no benchmark control is established for a particular dye, that plot will only display a black line that represents the spectrum of the unmixing control. The similarity index will display N/A.




- Similarity™ Matrix displays a similarity index matrix and a complexity index value.

Click View Similarity Index above the matrix to display the indexes for each dye. The Similarity™ Matrix will display the similarity index for each dye against itself and all the other dyes to be unmixed in the experiment.

The complexity index is a measure of how distinguishable a collection of spectral signatures are from each other when unmixed together. It calculates this by looking at the ratio of the

similarity index of the worst overlapping combination of signatures to the best overlapping combination of signatures.




10 Click Live Unmix.

Create New Unmixed Experiment


Live Unmix

11 The wizard closes and the experiment reappears. The reference group now has the unmixed icon to the left of the tube(s). Select an unmixed worksheet to view the unmixed data.



12 Select the sample tube you wish to acquire. The green arrow indicates the tube is selected. Click Start, then Record.

Use My Experiments to open experiments you ran if you wish to review the data or acquire more samples.



FCS files are stored in the Export folder by default, or the folder you set as the default. See [“Storage Preferences” on page 82](#) for information. FCS files for live unmixed data are saved as both raw data and unmixed data.

Analyzing Data Offline

To analyze data offline, you can click My Experiments, select the experiments you want to export, right-click and select Export. This will export the entire experiment as a ZIP file with all of the FCS files and worksheet templates contained inside. This experiment can be imported into other instances of SpectroFlo® CLC software, or unzipped to access the FCS files for analysis using other analysis software.

Post-Acquisition Unmixing

Samples can be acquired as raw data and then unmixed after acquisition is complete. This can be done through two methods:

- post-acquisition unmixing in the Acquisition module (see below)
- post-acquisition unmixing in the Extra Tools module (see [page 102](#))


Post-Acquisition Unmixing in the Acquisition Module

The unmixing wizard in the Acquisition module limits reference controls to those coming from the reference group in the experiment or reference controls run in QC & Setup.

To perform post-acquisition unmixing in the Acquisition module, perform the same workflow as live unmixing except the following:

- 1 Acquire all reference control tubes and sample tubes prior to selecting the Unmix button in the upper-left pane.
- 2 Examine the spectral plots by doing the following, if needed:
 - a. Move the polygon gate in the FSC vs SSC plot on the left to include the singlet population. Hold down Ctrl to move all the polygon gates at once.
 - b. Move the positive interval gate in the histogram to include the positively stained population. Move the negative interval gate to include the negative population.

- c. Move the interval gate on the spectrum plot on the right to select the channel that exhibits the brightest fluorescence intensity. This channel is the peak emission channel for the fluorescent tag.




- 3 Click Create New Unmixed Experiment.

A new experiment opens with a new unmixed worksheet.



Use My Experiments to open experiments you ran, if you wish to review the data or acquire more samples.



FCS files are stored in the Export folder by default, or the folder you set as the default. See ["Storage Preferences"](#) on page 82 for information. FCS files for post-acquisition unmixed data are saved as unmixed data only.

Post Acquisition Unmixing in the Extra Tools Module

When performing post-acquisition unmixing in the Extra Tools module, you can pick and choose which FCS files to unmix (for example, controls coming from different experiments, reference controls run during QC & Setup, or single-stained controls that were not run as part of the reference group).

FCS files can be designated into three categories:

- Single Stained
- Unstained
- Sample


■ **NOTE:** There must be at least one single-stained FCS file and one unstained FCS file in the file list. Otherwise, unmixing cannot be performed.

■ **NOTE:** The parameters need to match in order for this to work.


In addition, raw FCS files can also be conventionally compensated in this module through the Virtual Filters tab. This function can simulate the presence of filters and can compensate data using conventional compensation methods (see [“Virtual Filters” on page 106](#)).

To Unmix Raw Data Files:

- 1 Select Spectral Unmixing from the Extra Tools module.
- 2 Click Import to import raw FCS files for unmixing.



- 3 Select the files. Select multiple files using either the Shift or Ctrl key. Click **Open**.
- 4 Upon importing, a dialog box on how to assign sample types appears. Read the instructions and click OK.



- Once FCS files have been imported, select the sample type for each FCS file as Single Stained, Unstained, or Sample. The software will automatically designate the type based upon the file name. You can manually modify these if the automatic designation is incorrect.

Extra Tools Spectral Unmixing

Spectral Unmixing

Import Delete Clear All Show Plots

FCS File	Type	Fluorescent Tag	Label	☐ Universal Negative
Reference Group - FITC (Beads)	Single Stained			☐
Reference Group - PE (Beads)	Single Stained			☐
Reference Group - PerCP (Beads)	Single Stained			☐
Reference Group - Unstained (Beads)_L	Unstained			☐
Group_001 - MixedBeads	Sample			☐

- FCS files designated as single-stained will require a fluorescent tag designation to specify what reference spectrum will be provided for unmixing.

Extra Tools Spectral Unmixing

Spectral Unmixing

Import Delete Clear All Show Plots

FCS File	Type	Fluorescent Tag	Label	☐ Universal Negative
Reference Group - FITC (Beads)	Single Stained	FITC		☐
Reference Group - PE (Beads)	Single Stained	PE		☐
Reference Group - PerCP (Beads)	Single Stained	PerCP		☐
Reference Group - Unstained (Beads)_L	Unstained			☐
Group_001 - MixedBeads	Sample			☐

- Enter a label for each single-stained control and sample.

Extra Tools QC Beads

Spectral Unmixing

Import Delete Clear All Show Plots


FCS File	Type	Fluorescent Tag	Label	☐ Universal Negative
Reference Group - FITC (Beads)	Single Stained	FITC	CD4	☐
Reference Group - PE (Beads)	Single Stained	PE	CD8	☐
Reference Group - PerCP (Beads)	Single Stained	PerCP	CD3	☐
Reference Group - Unstained (Beads)_L	Unstained			☐
Group_001 - MixedBeads	Sample		4/8/3	☐

- Select Universal Negative for single-stained FCS files that do not contain a negative population, and the unstained control will be used for the negative population. In the bottom left of the screen, check whether Auto Fluorescence will be used as a fluorescent tag.

- 9 Click Show Plots to display the data in the FSC vs SSC plot, peak emission channel histogram, and spectrum plots.
- 10 The positive and negative populations need to be identified through the appropriate placement of the existing gates. Click OK to adjust the gates.
 - a. Move the polygon gate in the FSC vs SSC plot to include the singlet population.
 - b. Move the interval gate in the histogram labeled *Positive* to include the positively stained population. Move the interval gate in the histogram labeled *Negative* to include the negative population. Do not adjust the negative gate when using the Universal Negative.
 - c. Move the interval gate on the spectrum plot on the right to select the channel that exhibits the brightest fluorescence intensity. This channel is the peak emission channel for the fluorescent tag.



- 11 Click Unmix.
- 12 Select the directory to which the unmixed FCS files are exported or leave the default. Click OK.



These FCS files can then be imported to an experiment for analysis or analyzed using third-party software.


Adjusting Spillover

You can adjust the fluorescence spillover for a selected tube/well.

- 1 Ensure the green arrow is next to the tube/well and right click, then select Edit Properties in the group-tube hierarchy list.


The Tube/Well Properties window opens.

- 2 Select the Spillover tab.
- 3 Click the Enable Compensation checkbox to activate. The Adjust Spillover icon in the toolbar is now enabled.



- 4 Leaving the Tube/Well Properties window open, click the Adjust Spillover icon $\begin{bmatrix} 10 \\ 01 \end{bmatrix}$, then click and drag in the plot in the direction you want to adjust.

Changes are reflected in Tube/Well Properties window.




- 5 Click Save to save the changes. To discard the changes, close the Tube/Well Properties window and click No in the confirmation window, or click Discard changes in the tube/well properties window.

Virtual Filters

The Virtual Filters option in the Extra Tools module allows you to compensate raw FCS data using conventional compensation methods.


- 1 Click the Virtual Filters tab in the Extra Tools module.



- 2 Click Import to import raw FCS files for virtual filter analysis.

These FCS files can be single-stained reference controls, unstained controls, and/or sample files. However, you must include an unstained control FCS file.


- 3 Upon importing, a dialog box on how to assign sample types appears. Read the instructions and click OK.



- 4 Once FCS files have been imported, the sample type for each FCS file needs to be designated as Single Stained, Unstained, or Sample. The software will automatically designate the type based upon the file name. You can manually modify these if the automatic designation is incorrect.
- 5 FCS files designated as single stained require a fluorescent tag designation. Select the fluorescent tag for each single-stained sample.

If there is no negative population in the single-stained FCS file(s), select Universal Negative, and the unstained control will be used for the negative population.

The virtual filter is automatically assigned by the software based upon the fluorescent tag designation. (Optional) To increase the bandwidth of the virtual filter, use the channel pull-down menus to select the desired range. See the following table for wavelength ranges.




The following table shows the system's filter bandwidths.

Laser	Channel	Center Wavelength (nm)	Bandwidth (nm)	Wavelength Start (nm)	Wavelength End (nm)
Violet	V1	428	15	420	435
	V2	443	15	436	451
	V3	458	15	451	466
	V4	473	15	466	481
	V5	508	20	498	518
	V6	525	17	516	533
	V7	542	17	533	550
	V8	581	19	571	590
	V9	598	20	588	608
	V10	615	20	605	625
	V11	664	27	651	678
	V12	692	28	678	706
	V13	720	29	706	735
	V14	750	30	735	765
	V15	780	30	765	795
	V16	812	34	795	829
Blue	B1	508	20	498	518
	B2	525	17	516	533
	B3	542	17	533	550
	B4	581	19	571	590
	B5	598	20	588	608
	B6	615	20	605	625
	B7	661	17	653	670
	B8	679	18	670	688
	B9	697	19	688	707
	B10	717	20	707	727
	B11	738	21	728	749
	B12	760	23	749	772
	B13	783	23	772	795
	B14	812	34	795	829
Red	R1	661	17	653	670
	R2	679	18	670	688
	R3	697	19	688	707
	R4	717	20	707	727
	R5	738	21	728	749
	R6	760	23	749	772
	R7	783	23	772	795
	R8	812	34	795	829

6 (Optional) Select a label for the single-stained controls and samples.

7 Click Show Plots to display the plots.




FCS File	Type	Fluorescent Tag	Label	Universal Negative	Virtual Filter
Reference Group - FITC (Beads)	Single Stained	FITC		☒	Channel B2 to B2 - 518 - 539
Reference Group - PE (Beads)	Single Stained	PE		☒	Channel B4 to B4 - 560 - 583
Reference Group - PerCP (Beads)	Single Stained	PerCP		☒	Channel B8 to B8 - 669 - 687
Reference Group - Unstained (Beads)	Unstained			☐	
Group_001 - MixedBeads	Sample			☐	

The data is displayed in the FSC vs SSC plot and fluorescent tag histogram plot.


8 The positive and negative populations need to be identified through the appropriate placement of the gates. Click OK to adjust the gates.

- Move the polygon gate in the FSC vs SSC plot to include the singlet population. Hold down Ctrl to move all the polygon gates at once.
- Move the interval gate in the histogram labeled *Positive* to include the positively stained population. Move the interval gate in the histogram labeled *Negative* to include the negative population. Do not adjust the negative gate when using the Universal Negative.


The histogram x-axes are labeled with the fluorescent tag instead of the channel/detector.



- 9 Click Calculate Comp once gates have been set correctly. The conventionally compensated data is displayed. To view the data for a specific FCS files, select the file.



The spillover matrix is also calculated. Click Spillover Matrix to view the spillover values.




- 10 Click Export and select the location where you wish to export the conventionally compensated data. The files are exported to a folder named Compensated followed by the current date and time. Files can then be imported back into an experiment to analyze in SpectroFlo® CLC software.

Worklist

Create or open a worklist

- 1 In the Worklist module, click New in the toolbar to create a new worklist.
The existing worklists are displayed below the Worklist toolbar.
- 2 Select a desired Worklist and double-click, or click Open at the bottom to open the worklist.



Add tasks

After opening a worklist, click the **Add** button, or click the drop-down arrow in the Task column in the Worklist to add tasks.

Worklist_002						
Add Ref Controls Unmix Control Lot Options Approve All Export CSV Export PDF Print						
▼	No.	Sample ID	Task	Status	Signed By	
➔	1		Daily QC	▼		
➔	1.1		QC Tube (2002)			
▼	2	Sample_001	7C TBANK MT	▼	To Be Acquired	
	2.1		7C B			
▼	3	Sample_002	7C TBANK MT	▼	To Be Acquired	
	3.1		7C B			
▼	4		Fluidics Shutdown	▼		
	4.1		Fluidics Shutdown Tube 1			
	4.2		Fluidics Shutdown Tube 2			
	4.3		Fluidics Shutdown Tube 3			
	4.4		Fluidics Shutdown Tube 4			

- 1 Click **Add** to bring up the Add Task dialog box.
- 2 Select the Assay Task or Fluidics Task you want to add.

Add Task

Sample ID Start With:
 Start number of Sample ID:

▼ Assay Task	Carrier	Created By	Date Created	Description
➤ 7C TBANK MT	Manual Tube	Admin	June 22, 2021 - 14:48 PM	CytekAssaySetting
➤ 7C B MT	Manual Tube	Admin	June 22, 2021 - 14:44 PM	WB-LNW setting, Preview 10S, Med Flow rate
➤ 6C B TR	40 Tube Rack	Admin	June 16, 2021 - 00:01 AM	10S delay, Med Rate, 3000 Lym, 15 uL
➤ 6C B MT	Manual Tube	Admin	June 15, 2021 - 23:40 PM	Mid Rate Manual

▼ Fluidics Task	Description
Daily QC	Daily QC
SIT Flush	SIT Flush
Clean Flow Cell	Clean Flow Cell
Fluidics Shutdown	Fluidics Shutdown

Close
 Add

- 3 Enter the required number of samples in the quantity box in the lower right corner, and click the Add button.
- 4 When finished, click the Close button to close the dialog box.


Note that some fluidic tasks can only be added to specific carrier types. When adding a task, the software will automatically select the sample carrier required for the test and perform the default sample position layout on the carrier.

- 5 (Optional) Print the Tube/well location table. Click on Print and select Print Tube/Well Location Mapping Table.

Reference controls

You can choose to use the reference controls saved in the library or create reference controls in the worklist.

- 1 Click Ref Controls and the Create Reference Control dialog box will pop up.
The software will generate a reference control list according to the Assays added to the Worklist.
- **NOTE:** The Ref Controls icon will be grayed out until you add the assay you will be running.
- 2 Click the Add button to add additional reference controls.



Create Reference Controls

Reference Group Acquisition

Define Unstained Control(s) for Autofluorescence Extraction

Name: Unstained Control Type: Cells ☐ Define Additional Negative Control(s) for Spillover Calculation

Fluorescent Tag	Control Type	Label	Lot	Negative Control
FITC	Cells	CD3		
PE	Cells	CD16+56		
PE-Cy5	Cells	CD19		
PE-Cy7	Cells	CD4		
PE-Dazzle594	Cells	CD8		
PerCP-Cy5.5	Cells	CD45		

+ Add - Remove

Cancel Save Next

- 3 Click Acquisition or click Next to select worksheet, stopping criteria, and instrument setting for data collection.

If multiple assays in the Worklist use different user settings, all user setting plus CytekAssySetting will be listed in the drop-down menu.

CytekAssySetting will always be the default setting. The reference control data should be collected before the sample data.

✎ Create Reference Controls

Reference Group Acquisition

User Setting: 6CB ▾

CytekAssySetting (Cytek) ▾


Name	Sheet	Stopping Gate	Storage Gate	Events To Record	Stopping Volume (ul)	Stopping Criteria
Reference Group 6CB	Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
Unstained (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
CD45 FITC (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
CD56 PE (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
CD8 PE-Dazzle594 (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
PE-Dazzle594 (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
CD19 PE-Cy5 (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
PE-Cy5 (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
CD3 PerCP-Cy5.5 (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
CD4 PE-Cy7 (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
APC (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
CD19 APC (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
APC-Cy7 (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume
CD8 APC-Cy7 (Cells)	Default Raw Worksheet (Raw) ▾	All Events ▾	All Events ▾	5,000	3,000	☐ Count & Volume

Unmix

The Unmix process is the same as that in Experiment. See “[Unmixing Workflows](#)” on page 94 for details.

- 1 Click Unmix and select the corresponding reference control and whether to use autofluorescence extraction

2 Click Live Unmix.



Unmix Worklist

Select Controls

UNSTAINED CONTROLS

☒ Use Control from Library Negative (Cells)

☐ Use Control from Worklist

Name	Control Type

STAINED CONTROLS

☒ From Library	Fluorescent Tag	Control	Unstained	Generic
☒	FITC	CD3 FITC 20210301 (Cells)	▼	☒
☒	PE	CD56 PE (Cells)	▼	☒
☒	PE-Dazzle594	CD8 PE-Dazzle594 (Cells)	▼	☒
☒	PE-Cy5	CD19 PE-Cy5 (Cells)	▼	☒
☒	PerCP-Cy5.5	CD3 PerCP-Cy5.5 20210302-first (Cells)	▼	☒
☒	PE-Cy7	CD4 PE-Cy7 (Cells)	▼	☒
☒	APC	CD3 APC (Cells)	▼	☒


☒ Autofluorescence as a Fluorescent Tag

Cancel Live Unmix

Lot specific reference controls

You can select lot specific reference controls.

- 1 Click Control Lot to bring up the Reference control lot selection window.
- 2 If the same Fluorochrome has a different lot number, you can select a specific lot number in the drop-down menu. Not Specified Fluorochrome will use the default best fit reference control determined by the software.



Control Lot

6C B TR


	Sample ID	Task	FITC-CD3	PE-CD16+56	PE-Dazzle594-CD8	PE-Cy5-CD19	PerCP-Cy5.5-CD45	PE-Cy7-CD4
▼		WL-2-TR	▼	6C B(I) LNW ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼
▼ 1	Sample_001	6C B TR	6C BR LNW ▼	6C B(I) LNW ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼
1.1		6C B TR	6C BR LNW ▼	6C B(I) LNW ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼
▼ 2	Sample_002	6C B TR	6C B(I) LNW ▼	6C B(I) LNW ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼
2.1		6C B TR	6C B(I) LNW ▼	6C B(I) LNW ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼	Not Specified ▼

Cancel Save

Arrange sample locations in a carrier

Assays that contain carrier type cannot be changed in the Worklist, except assays that use Tube rack, which can be changed to manual in Worklist.

- 1 Click the Options button in the tool bar to open the sample carrier options dialog box.
 - You can change the sample carrier type (manual/Tube rack) if an assay uses the Tube rack.
 - You can customize sample tube/well locations on a carrier by drag and drop.



Data collection

- 1 Double-click an entry or a tube in the Worklist to make it the current task.
- 2 Click the Start button in the Acquisition Control panel and observe the scatter plot in the data display area on the right.
- 3 Adjust the FSC and/or SSC gain and Threshold in the instrument control pane to put the population of interest on scale. Adjust the gate in the FSC vs SSC plot to enclose the population.

- Click the Run button, and in the pop-up drop-down list, select the sample range to be collected to start the data collection. The software will prompt you to load the correct sample tube or Carrier type.

Worklist_003

Add
Ref Controls
Unmix
Control Lot
Options
Approve All
Export CSV
Export PDF
Print

No.	Sample ID	Task	Location	Carrier	Status
1	Sample_001	6C B TR	A1		To Be Acquired
1.1		6C B TR	A1	TubeRack_001	
2	Sample_002	6C B TR	A2		To Be Acquired
2.1		6C B TR	A2	TubeRack_001	
3	Sample_003	6C B TR	A3		To Be Acquired
3.1		6C B TR	A3	TubeRack_001	
4	Sample_004	6C B TR	A4		To Be Acquired
4.1		6C B TR	A4	TubeRack_001	

Acquisition Control

6C B TR: Not Ready - Warming up

Start
Run
Stop
Pause
Restart
Eject
SIT Flush

Flow Rate: Medium uL/Min: 0.00

Event Rate: 0 Abort Rate: 0

Threshold Count: 0

Time Elapsed: 00:00:00 (HH:MM:SS)

Events to Display: 5,000

Selection	Description
Run Current	Record data for current tube/well
Run Entry	Record data for current entry
Run From Entry	Start record data for the worklist from current entry
Run From Current	Start record data for the worklist from current tube/well
Run All	Start record data for the whole worklist from the first entry


Export Report

- Double-click an entry that has been recorded, and the report will be opened on the right side. The report enables you to preview the data for the following items: gating and unmixing accuracy, modify the report, add notes, approve and export reports.

Correct Unmix errors

If a worksheet or report was used for acquisition it is possible to modify compensation, if needed.

- 1 Right-click on a tube/well to and select Properties.
- 2 Go to the Spillover tab.
- 3 Adjust compensation and click Save.
- 4 The software will prompt you to select if you want to apply the compensation changes to other samples.



Modify Report

The components on a Report cannot be added or removed but the gates can be adjusted according to the data.

- 1 After the adjustment is completed for one entry, you can right-click on the entry (sample) to apply the adjustment to other samples in the worklist.

The screenshot displays the 'Worklist' window with a table of sample entries. A right-click context menu is open over the entry 'Sample_001' (Task: 6CB). The menu options are: Insert Task(s), Delete, Run Selected, Run All, Paste Spillover, Apply Report (highlighted), Approve, and Unapprove. A tooltip for 'Apply Report' states: 'Apply changes in "Report_001" to other tasks'. The background shows the 'Acquisition Control' panel with 'S1: Ready' and a flow rate of 'Low'. A scatter plot titled 'All Events' is visible on the right, showing a gate labeled 'P1'.

No.	Sample ID	Task	Status	Sign
1.4		CD8 PE-Dazzle594 (Cells)		
1.5		CD19 PE-Cy5 (Cells)		
1.6		CD45 PerCP-Cy5.5 (Cells)		
1.7		CD4 PE-Cy7 (Cells)		
2	Sample_001	6CB	To Be Approved	
2.1		S1		
3	Sample_002	6CB		
3.1		S1		

Audit trail

The software will automatically record the key operations and manual changes in the audit trail log.

- 1 Click the Audit Trail Logs button to open the Audit Trail Logs dialog box.

The screenshot shows the top navigation bar with the following buttons: QC & Setup, Acquisition, Worklist (active), Extra Tools, Library, Preferences, Users, Help, and Sign Out. The status bar at the bottom indicates 'TUBES 6C B TR(6C B TR)' and 'REPORTS 6C B TR'. On the right side of the status bar, there are buttons for 'Approve' and 'Audit Trail Logs'.

- 2 In the audit trail log dialog box, you can enter the reasons in batches, or one by one.

Audit Trail Log

Export All Audit Trail Logs

No.	Action	Reason	Operator	Date Created
2	6 CBR Unapproved by Administrator at June 11, 2021 - 13:15 PM.			
	Approved by Administrator at June 09, 2021 - 09:33 AM. Signature Comments: test		Administrator	June 09, 2021 - 09:33 AM
	Unapproved by Administrator at June 09, 2021 - 09:32 AM. Comments:		Administrator	June 09, 2021 - 09:32 AM
	Approved by Administrator at June 09, 2021 - 09:29 AM. Signature Comments: test		Administrator	June 09, 2021 - 09:29 AM
	Modify Gate in Report: Report_001.		Administrator	June 09, 2021 - 09:29 AM
	Modify Gate in Report: Report_001.	Modify	Administrator	June 09, 2021 - 09:29 AM
	Modify Gate in Report: Report_001.		Administrator	June 09, 2021 - 09:29 AM
	Modify Gate in Report: Report_001.		Administrator	June 09, 2021 - 09:29 AM
	Modify Gate in Report: Report_001.		Administrator	June 09, 2021 - 09:29 AM
	Modify Gate in Report: Report_001.		Administrator	June 09, 2021 - 09:29 AM
	Modify Gate in Report: Report_001.		Administrator	June 09, 2021 - 09:29 AM

Reason

Modify

Close

Approve Report

- 1 After report adjustment and review are completed, click the Approve button to approve the report.
- 2 The approval process will prompt for the user name and password of the approver, and allow comments.
- 3 After completion, the status of the test task will change to Approved.

- Click **Approve All** for batch approval after confirming that all reports are correct.

Worklist Worklist Entries						
Worklist_003						
+ Add Ref Controls Unmix Control Lot Options Approve All Export CSV Export PDF Print						
No.	Sample ID	Task	Location	Carrier	Status	
1	Sample_001	6C B TR	A1		To Be Approved	
1.1		6C B TR	A1	TubeRack_001		
2	Sample_002	6C B TR	A2		To Be Approved	
2.1		6C B TR	A2	TubeRack_001		
3	Sample_003	6C B TR	A3		To Be Approved	
3.1		6C B TR	A3	TubeRack_001		
4	Sample_004	6C B TR	A4		To Be Approved	
4.1		6C B TR	A4	TubeRack_001		

- If you want to re-record or make changes in an approved report, click the **Unapprove** button.
- The entry status will change from approved to pending approval.
- You can re-record, or adjust the gates, etc., and re-approve the entry.

Print Report

After a report is approved, the PDF file of the report and CSV file of the data in statistical table will be automatically exported to C:\CytekbioExport_CLC\WorklistReportResults. You also can manually export PDF file and CSV file to a selected directory.

Worklist_003						
+ Add Ref Controls Unmix Control Lot Options Approve All Export CSV Export PDF Print						

- To print Reports, click **Print**.
- Select the test report to be printed and the audit trail log to be included.

3 Follow the rest of the prompts to print.

 **Print**


☒	No.	Task	Status	☒ Include Audit Trail Logs
☒	1	 6C B TR	Approved	☒
		 6C B TR		
☒	2	 6C B TR	Approved	☒
		 6C B TR		
☒	3	 6C B TR	Approved	☒
		 6C B TR		
☒	4	 6C B TR	Approved	☒
		 6C B TR		

 Cancel  **Print Preview**

Loader

Automated Sample Loader Components


The loader power switch is located on the back of the loader. To power up the loader, connect the power supply and turn the power switch to the ON position.




When the loader is powered on and successfully connected to the software, the status indicator in the lower-right corner of the screen shows a green checkmark. A red X appears if the loader is not powered on or not connected to the software.

✓ Sheath ✓ Waste ✓ Cytometer **✓ Loader**



- The right panel includes the main controls for the loader. There is one USB Port, a power supply connector, and a power switch.



- The loader has 2 covers that can be opened to put the sample carrier onto the stage platform.
- The fan is located on the right side and is used for ventilation of the electronics. It comes on when the temp is above 27deg C.
- The software detects the position of the handle and will set the cytometer to function in the corresponding mode.




- The Manual Handle is located on the right side of the loader. The handle is a mechanical lever used to make the instrument work in Plate Mode, Tube Mode, or Tube Rack Mode. The handle can be moved with either power on or power off.



Loader Acquisition Controls and Settings

The Acquisition Control pane allows you to start, stop, and pause acquisition, record data, and restart acquisition counters. To show, hide, or undock (float) this pane from the experiment panel, use the dock/undock and hide icons in the top-right corner.



The following table describes the controls in the Acquisition Control pane.

No.	Control	Description
1	Start/Record/Stop/Load/Eject/Restart	It will depend on whether running a plate or running a tube rack. Running a plate: the lever is forward. Running a tube rack: the lever is to the back. Select Start after selecting Load plate to start acquisition. Select Record to record data. Record can also start acquisition. Select Stop to stop acquisition. Select Load or Eject to load or eject the plate. Select Restart to restart the acquisition counters. All events and results displayed are refreshed. Stop and Restart are enabled once Start is selected. Stop and Pause are enabled once Record is selected.
2	SIT Flush	Select to perform a SIT Flush
3	Flow Rate	Select Low (15 µL/min), Medium (30 µL/min), or High (100 µL/min). The exact flow rate is displayed.
4	Event Rate, Abort Rate, Threshold Count, Time Elapsed	Displays the real-time counts during acquisition.
5	Events to Display	Enter the number of events to display during acquisition.

The following table describes the ASL settings.

Setting	Description
Selected Settings	Three loader settings are available—default, high throughput, and low carryover. You can create your own custom settings.
Acquisition order	Select the order that you want the plate to run. Wells are acquired by: • row from left to right (A1-A12, B1-B12, etc.) • column from top to bottom (1A-1H, 2A-2H, etc.) • row from left to right, then right to left (A1-A12, B12-B1, C1-C12, etc.) • column from top to bottom, then bottom to top (1A-1H, 2H-2A, etc.)
Shake Time	Select the time (in seconds) to shake the plate/tube rack. You can also disable shake time.
Shake Speed	Select the speed of the orbital shaker (in RPM).
Shake Interval Mode	Select whether you would like to shake every N wells/tubes or after a specified period of time. You can also disable shaking.
Shake Every N Wells, or Shake Interval	Select how often (number of wells/tubes or time in seconds) to shake the plate/tube rack.
Premix Time	Select the time (in seconds) to shake the plate/tube rack before acquisition of the first tube/well.

Setting	Description
SIT Flush Times	A SIT flush is performed over the wash station after each acquisition. Choose single flush, double flush, or disabled if you do not wish to perform a SIT flush.
Sample Recovery	Allows any remaining sample that is left in the SIT after acquisition is complete to be deposited back into the well.
Record Data Delay Time	Select the time in seconds you wish to preview data from a well/tube before recording begins once you click Record.

Loader Daily Workflow

The following section describes the recommended procedures for operating the Automated Sample Loader from start up to shut down.

Startup

- 1 Power on the cytometer, loader, and computer (if not already powered up).
- 2 Launch the SpectroFlo® CLC software and login with the appropriate username.
- 3 Once logged into the software, the cytometer initialization procedure will begin. Load a tube of DI water and wait for the instrument and the loader self-checking procedure to finish. The Cytometer and loader connection status indicator will turn green when self-checking is complete.
- 4 Ensure that the status indicator checkmark for “Cytometer” and “Loader” are both green. This indicates that the workstation, instrument, and loader are communicating with one another.



- 5 The loader and cytometer are now ready for use.

Set Up Experiment

Setting up the experiment in SpectroFlo® CLC includes:

- 1 Providing a name and description (optional) for the experiment.
- 2 Specifying the fluorophores used in the experiment.
- 3 Defining the reference group as either acquired in the experiment or referred from the stored reference control library.
- 4 Defining the acquisition criteria for reference group and sample group. Acquisition criteria is based on events, volume or time.
- 5 (optional) If the reference controls are to be acquired in the experiment, it CANNOT be live unmixed during acquisition. Acquisition must proceed with a raw worksheet. Reference controls that are referred to from the stored reference control library do not need to be acquired in the experiment.
- 6 Run the sample, adjust instrument setting, and save User Settings in the software.


Acquisition

The Acquisition workspace provides the user with the necessary elements for data collection. Flow cytometer data can be acquired from experiments. Experiments can be created either through the new experiment wizard or using experiment templates. SpectroFlo® CLC software includes many features that ensure the proper acquisition of flow cytometry data. Data is recorded into the FCS 3.0 format.

The Acquisition workspace provides an interface where the user can layout an experiment worksheet with different elements such as 2-D plots, histograms, and statistics displays. Instrument gain settings are also adjustable in the interface. Properly adjust instrument gain settings between samples, and select whether to preview or record data with the appropriate acquisition criteria at selected sample flow rates.

Creating a new Loader experiment with Reference Group:

- 1 Select Edit in the Acquisition menu to edit the Plate experiment.



- 2 The new experiment wizard opens. In this window, the user can specify a name for the experiment and/or type in a description.

 Edit Experiment

Fluorescent Tags Groups Markers Keywords Acquisition Loader Settings

Name

Experiment_001

Description (optional)

Default Experiment

- 3 The user must then select the fluorophores that are present in the experiment from the Library list pane on the left. The user needs to select all individual fluorophores present in the experiment as this will determine which reference controls are to be used during spectral unmixing. Additional fluorophores that are not present in the sample can also be selected from this list.

Library

▼ Blue Laser

BB515

Alexa Fluor 488

FITC

Alexa Fluor 532

PE

PE-Dazzle594

PE-eFluor 610

PE-CF594

PE-Texas Red

PE-Alexa Fluor 610

PE-Cy5

PerCP

PE-Cy5.5

PerCP-Cy5.5

BB700

PerCP-eFluor 710

PE-Alexa Fluor 700

PE-Cy7


► Red Laser

► Violet Laser


► Fluorescent Proteins

► Viability


■ **NOTE:** You can search for fluorophores in the upper right corner of the Library pane.




- 4 Once all fluorophores have been chosen from the fluorophore library, confirm in the selection pane. Individual fluorophores can be removed in this pane. Select **Clear All** to remove all fluorophores.




- 5 Select the wells on the plate to designate them for acquisition. Configure them as either a sample, reference, or clean well.



- 6 Using the Acquisition tab, configure the parameters of the experiment. These include display worksheet, stopping and storage gates, and acquisition criteria.



Using the **Loader Settings** tab, configure the operation parameters of the loader. These include shake configuration, SIT Flush Times, Sample Recovery, and Record Data Delay Time.



Loader Shutdown



Always adhere to proper handling procedures for samples and reagents. Wear standard PPE laboratory attire such as protective gloves, eyewear, and lab coat during this procedure.

Fluidics shutdown with loader tube rack:

- 1 In the Cytometer Wizard, select Fluidics Shutdown with Tube rack.
- 2 Load the 10% Bleach tube into the A1 position (Tube position can be changed by dragging on plate view).
- 3 Load the deionized water tube into B1 position (Tube position can be changed by dragging on plate view).
- 4 Load the 30% Contrad tube into C1 position (Tube position can be changed by dragging on plate view).
- 5 Load the deionized water tube into D1 position (Tube position can be changed by dragging on plate view).

- 6 Leave the SIT submerged underneath deionized water of the last tube.
- 7 Exit SpectroFlo® CLC software.
- 8 Turn off the cytometer and computer.
- 9 Turn off the loader.

Maintenance

Maintenance Schedule



Any instrument surface in contact with biological specimens can transmit potentially fatal disease. Use universal precautions when cleaning the instrument or replacing parts. Wear suitable laboratory attire such as protective gloves, eyewear, and lab coat.

The 10% bleach solution used throughout the maintenance procedures is prepared by adding 1 part household bleach to 9 parts DI water.

Routine maintenance of the Cytex Northern Lights™-CLC cytometer includes periodic replacement of parts. For part numbers, see [“Supplies and Replacement Parts” on page 161](#).

Scheduled Maintenance

The following table describes the scheduled maintenance procedures for your cytometer.

Maintenance Procedure	Description	Frequency
Replace sheath filter	Ensures debris-free sheath fluid	Every 6 months, or as needed
Long Clean	Cleans the fluidic lines with 10% bleach solution	Once a month and prior to service calls

Unscheduled Maintenance


The following table describes the unscheduled maintenance procedures for your cytometer.

Maintenance Procedure	Description	Frequency
SIT Flush	Backflushes the SIT	As needed, if SIT clogs or after running sticky dyes
Purge Filter	Removes bubbles from the sheath filter	If bubbles are present in the sheath filter, or if the plenum or sheath tank runs dry
Clean Flow Cell	Runs 10% bleach solution followed by DI water through the flow cell	As needed, or after running sticky dyes
Clean external surfaces	Keeps surfaces free from salt buildup	As needed

SIT Flush

A sample line backflush is performed whenever a tube is removed from the SIP after sample acquisition. If the sample line exhibits signs of carryover or becomes clogged after completing an experiment with a sticky dye such as propidium iodide, acridine orange, or thiazole orange, the sample line should be manually backflushed.

- 1 In the Cytometer tab, from either the QC & Setup or Acquisition module, select SIT Flush.



- 2 If carryover or a clog persists, place a tube of 10% bleach on the SIP and acquire at High flow rate for 5 minutes. Afterwards, acquire a tube of DI water at High flow rate for 5 minutes.


■ **NOTE:** If running large quantities of beads or large cells, we recommend running a tube of 10% bleach followed by a tube of DI water, each for 5 minutes, between experiments.

Purge Filter



Perform this procedure if air bubbles are visible in the sheath filter, or if the plenum or sheath tank have run dry and air is present in the fluidics system.

- 1 In the Cytometer tab, from either the QC & Setup or Acquisition module, select Purge Filter.



The vent valve connected to the sheath filter will open releasing any air bubbles trapped inside the sheath filter.


- 2 Repeat the Purge Filter fluidic mode until there are no visible bubbles inside the sheath filter.

Clean Flow Cell

Clean the flow cell after completing an experiment with a sticky dye such as propidium iodide, acridine orange, or thiazole orange. Cleaning the flow cell is also recommended after acquiring large quantities of highly concentrated bead solutions or if you suspect a clog.

If running sticky samples, perform this procedure using a 50% Contrad 70 solution instead of 10% bleach between experiments.

- 1 In the Cytometer tab, from either the QC & Setup or Acquisition module, select Clean Flow Cell.



- 2 Follow the instructions that appear. Load a tube containing 3 mL of 10% bleach solution on the SIP and click Continue.
- 3 When prompted, load a tube containing 3 mL of DI water on the SIP and click Continue.
- 4 Click Done when the procedure is complete.


Long Clean

Decontaminate the fluidics system monthly by running the Long Clean fluidics mode. Run the Long Clean just prior to service calls, or if you run high volumes of unwashed samples or samples stained with propidium iodide, acridine orange, or thiazole orange.



Do not run bleach or detergent through the sheath filter. It is difficult to remove cleaning solutions from the sheath filter.

- 1 In the Cytometer tab, from either the QC & Setup or Acquisition module, select Long Clean.



- 2 Follow the instructions that appear. Prepare the appropriate cleaning tubes and fluidics tanks.
- 3 Empty the waste tank. Replace the sheath filter with the sheath filter bypass (long clean tubing) assembly.



- 4 Detach the sheath tank and replace it with a tank containing a 10% bleach solution.
- 5 Load a tube containing 3 mL of a 10% bleach solution on the SIP.
- 6 Proceed with the Long Clean in the software.
- 7 Once the bleach cleaning cycle is complete, reattach the sheath tank.
- 8 Remove the tube of 10% bleach from the SIP and replace with a tube of 3 mL of DI water.
- 9 Proceed with the Long Clean in the software.
- 10 When prompted, remove the long clean tubing assembly and re-install the sheath filter.

Fluidics Shutdown

Perform fluidic shutdown at the end of each day that you use the instrument. The shutdown procedure thoroughly cleans the fluidics system.

- 1 In the Cytometer menu from either the QC & Setup or Acquisition module, select Fluidics Shutdown.
- 2 Load a tube containing 3 mL of 10% bleach on the SIP and click Continue.
A 10% bleach solution is prepared by adding 1 part household bleach to 9 parts DI water.
- 3 Load a tube containing 3 mL of DI water and click Continue.
- 4 Load a tube containing 3 mL of 50% Contrad 70 and click Continue.
- 5 Load a tube containing 3 mL of DI water and click Continue.
- 6 Allow the shutdown procedure to complete, then click Done and turn off the cytometer. Make sure the SIT is submerged in the DI water at the end of the procedure.

The following day when you turn on the system, the startup procedure begins, using the tube of DI water on the SIP.

Cleaning the External Surfaces

Periodically check for saline residue.

- 1 Dampen a cloth with a mild cleaning solution and wipe the surfaces of the instrument.
- 2 Dampen a cloth with DI water and wipe the surfaces again to remove residual cleaning solution.
- 3 Dry the surfaces with a clean, dry cloth.

Inspecting the Fluidics Lines



Check the cytometer periodically for fluid leaks. If any evidence of a leak is detected, contact Cytex Technical Support immediately. Do not attempt to repair the instrument.

- 1 Visually inspect for fluid leaks by looking for small pools of liquid near any of the quick-connects.
- 2 Visually inspect for dried residue or slight discoloration in the spaces surrounding the cytometer.

Replacing the Sheath Filter


The sheath filter traps debris and air bubbles before they reach the flow cell. Replace the filter assembly every 6 months, or when you see increased sample flow rates or debris in an FSC vs SSC plot.



Wear appropriate safety attire such as protective gloves, eyewear, and lab coat while performing this procedure.

- 1 Turn off the cytometer.
- 2 Open the front cytometer panel.

- 3 Press the two fluidics line quick-connects and the vent line quick-connect to the right of the sheath filter.



- 4 Discard the sheath filter according to standard laboratory protocol and local regulations.
- 5 Install a new sheath filter with the arrow pointing up.
- 6 Restart the cytometer and run the Purge Filter fluidic mode to remove air bubbles (see [“Purge Filter” on page 136](#)). Repeat this step until all air bubbles are purged from the filter.
- 7 Close the front panel.

Replacing the SIT


Replace the SIT if the tubing is clogged even after repeatedly cleaning and flushing the SIT.

- 1 Ensure the SIT is extended in a tube of DI water and the cytometer is turned off.

If the cytometer was shut down properly using the Fluidics Shutdown procedure, the SIT will already be extended and left in water.


■ **NOTE:** If the SIT is not extended, turn on the cytometer and run Fluidics Shutdown. Then turn off the cytometer.

- 2 Obtain a SIT tubing assembly.



- 3 Open the SIT door. Identify the three components shown in the following figure.

Open the front panel, then open the SIT door. The SIT door is located above the SIP. See [“Front of Cytometer”](#) on page 16.



- 4 Twist off and carefully remove the black plastic nut from the bottom of the flow meter.



- 5 Follow the tubing from the black nut down to the beige plastic nut. Twist off the beige nut and gently pull it and tubing out from the SIP.




- 6 Discard the SIT tubing assembly.

- 7 Insert the new sample tubing through the hole.



- 8 Slide the beige nut down into the hole and turn the nut to secure, but do not tighten.




- 9 Mark the fully screwed in position at the base and on the nut, as shown. Then unscrew the nut $1\frac{1}{2}$ revolutions from the marked position.




Mark position of screwed in nut.

- 10 Verify that the tubing contacts the bottom of the tube on the SIP.



- 11 Secure the black nut to the bottom of the flow meter so it is firmly attached. Make sure the sample line bends away from the SIT assembly when the SIT stage rises.



Secure black nut.

Ensure sample line bends away from SIT assembly.

- 12 Close the SIT door.
- 13 Turn on the system. During initialization the system automatically calibrates the SIT depth in the tube loaded on the SIP.

Loader Maintenance




Any instrument surface in contact with biological specimens can transmit potentially fatal disease. Use universal precautions when cleaning the instrument or replacing parts. Wear suitable laboratory attire such as protective gloves, eyewear, and lab coat.

Calibrating the Loader Stage and Plate

Click the Reset button to recalibrate the stage and plate if either of the following occurs:

- If the stage itself reaches the limit of its mechanical travel range.

- If any component that moves along with the stage reaches the limit of its mechanical travel range.



- 1 Pull the handle to the front. Observe the tube guide position and adjust the tube guide maximum travel to a suitable position.
- 2 Click Alignment Point A or Alignment Point B. Align the tube guide with the cross-shaped marker of the stage calibration fixture in the left/right and back/forth positions.
- 3 Place a standard well plate or deep well plate onto the shaker stage. Calibrate tube guide position and A1 position.
- 4 Pull the handle back and place the tube rack onto the shaker stage. Calibrate A1 position.


Cleaning the External Surfaces of the Loader

Periodically check for saline residue.

- 1 Dampen a cloth with cleaning solution and wipe dried sheath fluid from surfaces.
- 2 Wet a clean cloth with DI water and wipe surfaces to prevent bleach corrosion.
- 3 Wipe surfaces with a clean, dry cloth.

■ **NOTE:** High concentrations of sodium hypochlorite (bleach) and other cleaning agents can damage the instrument.

Loader Connections



- USB2.0: Requires a standard USB2.0 type-A transfer USB2.0 type-A cable (approx. 1 meter length) for connection to computer.
- Power: Requires an IEC320-C8 power cable to connect with AC-DC power adapter.

Uninstalling the Loader Power Supply

To uninstall the power supply:

- 1 Unplug the AC connection.
- 2 Unplug the loader power supply connection.



Loader Lifting and Carrying Instructions

- Keep all equipment off the ground on a firm surface such as a table.
- Never place the loader on the ground when in use.
- Use a cart or other stable support system to carry the loader.

Troubleshooting

This section provides tips to help you identify and resolve issues that might occur on your flow cytometer. If additional assistance is required, contact Cytex Biosciences. Please have the following information available: serial number, error messages, and details of recent performance.

For instrument support within the US, call 1-877-92-CYTEK. Visit our website, www.cytexbio.com, for up-to-date contact information.

- “General Troubleshooting” on page 148
- “Loader Troubleshooting” on page 150

General Troubleshooting

Observation	Possible Causes	Recommended Solutions
Daily QC does not complete	Wrong QC bead sample	Ensure you are running SpectroFlo® QC Beads.
	Bead sample not properly mixed	Mix the bead sample.
	Bead sample too dilute	Concentrate the bead sample or prepare a fresh bead sample.
	Air bubble in sample line	Run a SIT Flush.
Daily QC failed	Air bubble in fluidics	Run a Purge Filter.
	Dirty flow cell	Run a Clean Flow Cell. If the problem persists, run a Clean Flow Cell using 25%–50% Contrad 70, followed by DI water.
	Questionable sample prep	Verify the sample prep technique.
	Air in sheath filter	Run a Purge Filter.
	Sample not diluted in same fluid as sheath	Dilute the sample in the same fluid as the sheath solution.
Air in sheath filter	Cytometer was not in use for a prolonged period	Run a Purge Filter. Check that all sheath connectors are securely attached. Check for leaks or cracks in the sheath plenum. Replace, if needed.
	Empty sheath tank	Fill the sheath tank. Run a Purge Filter.
No events displayed (flow rate lower than expected)	No sample in tube	Add sample or install a new sample tube.
	Sample not properly mixed	Mix the sample to suspend cells/particles.
	Clogged SIT	Run a SIT Flush. Then run a Clean Flow Cell with 10% bleach, followed by a Clean Flow Cell with DI water. If the clog persists, replace the sample line.

Observation	Possible Causes	Recommended Solutions
No events displayed (flow rate normal)	Insufficient gain for threshold parameter	Increase the gain for the threshold parameter.
	Threshold too high	Lower the threshold.
	Laser delay not correct	Ensure the laser delay values match those from the latest Daily QC run. See “Instrument Control” on page 50 for the laser delay location. If the values do not match, rerun Daily QC.
	Threshold set to incorrect parameter	Set the threshold to the appropriate parameter for the application (usually FSC).
	Gated plot with no data in gate	Delete or move the gate.
Low sample event rate	Threshold too high	Lower the threshold.
	Insufficient gain for threshold parameter	Increase the gain for the threshold parameter.
	Sample not properly mixed	Mix the sample to suspend cells/particles.
	Sample too dilute	Concentrate the sample. Set the flow rate to Medium or High.
	Clogged SIT	Run a SIT Flush. Then run a Clean Flow Cell with 10% bleach, followed by a Clean Flow Cell with DI water. If the clog persists, replace the sample line.
Erratic event rate	Partially blocked SIT	Run a SIT Flush. Then run a Clean Flow Cell with 10% bleach, followed by a Clean Flow Cell with DI water.
	Clumpy sample	Vortex, filter, or disaggregate the sample.
Data in scatter parameters appear distorted	Air bubble in flow cell	Run a SIT Flush.
	Air in sheath filter	Run a Purge Filter.
	Dirty flow cell	Run a Clean Flow Cell.
	Poor sample health	Check the viability of the cells.
	Hypertonic buffers	Check the pH of the buffers and fixative.
	Incorrect instrument settings	Optimize the instrument settings.

Observation	Possible Causes	Recommended Solutions
High CVs	Air bubble in fluidics	Run a SIT Flush and a Purge Filter.
	Sample flow rate set to High	Set the sample flow rate to Low or Medium.
	Dirty flow cell	Run a Clean Flow Cell. If the problem persists, run a Clean Flow Cell using 25%–50% Contrad 70, followed by DI water.
	Questionable sample prep	Verify the sample prep technique.
	Air in sheath filter	Run a Purge Filter.
	Sample not diluted in same fluid as sheath	Dilute the sample in the same fluid as the sheath solution.

Loader Troubleshooting

Observation	Possible Causes	Recommended Solution
Stage or motor not moving	USB connection	Check that the USB cable and USB connector are not loose.
	Power connection	- Check that the Adapter power, AC-DC power cable, and DC-DC power cable have good connections. - Use a meter to check if the DC output of the adapter power is 24VDC.
	Mechanical breakdown	With power switch OFF, manually push and pull the stage to check if the motion is smooth.
Stage is not aligned	Stage calibration is not correct.	- Re-calibrate the instrument by following the calibration steps in software. - Check if there is significant skew on the stage. - Check if the software configuration for Plate/Deep Well/Tube Rack is correct.

Glossary

ASL	Automated Sample Loader mixes samples and delivers tube racks and plates to the cytometer for acquisition.
APD	Avalanche photodiode (APD) is a highly sensitive semiconductor electronic device for measuring light intensity. APDs convert light to electricity and can be thought of as photodetectors.
auto-fluorescence	The inherent fluorescence arising primarily from cell structures such as mitochondria and lysosomes. Auto-fluorescence can hinder detection of dim fluorescent signals.
compensation	The process by which spillover fluorescence from secondary parameters is accounted for so that fluorescence values for a parameter represent only the fluorescence of the primary fluorophore.
data file	A collection of measured values from a single sample combined with text describing the data that has been stored as a flow cytometry standard (.fcs) file to disk.
deconvolve	An algorithm-based process used to reverse the effects of convolution (or overlapping) on recorded data.
detector	A device that responds to a specific stimulus. Photodiodes and APDs tubes are two types of detectors in cytometers. They convert light signals into electronic signals.
dot plot	A graphical representation of two-parameter data. Each axis of a plot displays values of one parameter.
electronic noise	Random fluctuation in electronic signals, a characteristic of all electronic circuits.
event rate	The rate at which cells or particles are acquired.

FCS	Flow cytometry standard, a standard format for flow cytometry data files.
filter	An optical device that blocks the passage of part of the incident light, allowing the rest to pass virtually unchanged.
flow cell	The flow cell enables hydrodynamic focusing of the sample so that the individual cells or particles of interest can be interrogated by the laser(s) sequentially.
flow cytometry	A technology that simultaneously measures and analyzes multiple characteristics of single cells or particles as they pass through a laser beam.
flow rate	The amount of fluid passing through a point per unit of time. The Cytex Northern Lights™-CLC cytometer's flow rate is measured by the flow meter.
fluorescence	The emission of light of longer wavelengths that occurs when a substance absorbs light of shorter wavelengths.
fluorophore	A fluorescent dye. A molecule capable of absorbing light energy, then emitting light at a longer wavelength (fluorescence) as it releases this energy.
gain	Amplification of a signal. Increasing gain results in a larger output signal for a given input signal.
gate	A numerical or graphical boundary (region) that defines a subset of data. Gates can be single- or multi-dimensional.
laser	Light Amplification by Stimulated Emission of Radiation. A light source that is highly directional, monochromatic, coherent, and bright. The emitted light is in one or more narrow spectral bands, and with most lasers, is concentrated in an intense, narrow beam.
laser delay	Amount of time between signals from different laser intercepts.
photodiode	A device for measuring light intensity. A photodiode generates an output current proportional to the incident light intensity.
rCV	Percent robust coefficient of variation. The robust CV is a measure of dispersion of a probability distribution or a frequency distribution.

rSD	Robust standard deviation. The robust SD is based on the deviation of individual data points to the median of the population.
reference	Spectral profile of a fluorescent tag in all detectors for all lasers.
resolution	A measure of a cytometer's ability to distinguish between two populations with differing fluorescence or light scatter intensities.
sensitivity	A measure of a cytometer's ability to distinguish particles from background noise. It is often expressed in terms of a minimum number of fluorescent molecules per particle required to distinguish a stained particle from an unstained particle. Sensitivity depends on the instrument, the dye, and the preparation method.
SIP	Sample injection port. The area of the cytometer where the sample is placed.
SIT	Sample injection tube. The probe that pulls sample from the sample tube to the flow cell.
spectral overlap	The phenomenon of different fluorophores emitting light within the same detection range. In multi-color experiments, compensation must be performed to correct for spectral overlap.
spectral unmixing	The mathematical method of deconvolving signals from multiple fluorophores/dyes to differentiate them.
spillover	Emitted light from a fluorophore entering the detector of another fluorophore.
voltage	Measure of electric potential. The voltage applied to a detector affects its amplification gain.

Specifications

Description	Specification
Workstation	
Operating system	Microsoft® Windows® 10 Professional x64 or Windows® 11 Professional
Processor	Intel® Core™ i5-12500 or equivalent
Memory	16GB RAM, 4800 DDR5 or higher
Storage	1TB SSD
Network Adapter	802.11ac 2.4/5 GHz wireless adapter
Monitor	32-in UHD 4K Monitor
Video Card	Integrated, additional I/O port card DP required
Installation Requirements	
Recommended workspace	152.4 x 61 x 132 cm (60 x 24 x 52 in)
Power	110/240 V~, 50/60 Hz, 3.5 A/1.5 A
Fuse rating	250 VAC, 5 A, size 5 x 20-mm
Heat dissipation	500 W with all solid-state lasers
Temperature	15°C–28°C
Environmental conditions	Indoor use
Altitude	Up to 2,000 m (6,562 ft)
Pollution degree	2
Humidity	20% to 85% relative non-condensing
Air filtering	No excessive dust or smoke
Lighting	No special requirements
General Parameters	

Description	Specification
Instrument Composition	The device is made up of the main unit (fluidics subsystem, optical subsystem, electronic subsystem), fluid reservoirs, software and the optional ASL.
Instrument Function	- Data collection, analysis, storage, and output - Analytic capability for FSC and SSC - Analytic capability for up to 38 fluorescence channels - Capability of measuring the acquired sample volume - Automatic sample loading with tube-rack or plate using ASL
Dimensions (W x D x H)	Main Unit: 54 x 52 x 52 cm (21.3 x 20.5 x 20.5 in) ASL: 23 x 47 x 28 cm (9.1 x 18.5 x 11 in)
Weight	Main Unit: 61 – 71 kg (134.5 – 156.5 lb.) ASL: 15 kg (33.1 lb.)
Software	
Name & Version	SpectroFlo® CLC v1.0.4 or above
Run Environment	Microsoft Windows® 10 or Windows® 11 Professional (64-bit)
Software Function	- User management: 3-tier user accounts (administrator, lab manager and lab technician) - Daily QC: automated routine performance tracking and set up of NL-CLC - Data acquisition with live unmixing using experiment and worklist - Instrument control: acquisition settings (threshold, events count, parameter selection), ASL settings (shake speed, shake time, shake interval, SIT flushes), acquisition modes (default, low carryover, high throughput) - Data analysis with audit trail - Data export: experiments, worklists, raw and unmixed FCS files, statistics, templates
Data Acquisition Rate	The system is capable of acquiring up to 35,000 events/s
Normal Operating Conditions	
Temperature Range	Operation: +15°C to +28°C Transportation & Storage: -20°C to +55°C
Relative Humidity	Operation: 20%-85% non-condensing Transportation & Storage: 10%-90%
Power Supply	Main Unit: 110-240 V AC, 50/60 Hz, 3.5-1.5 A ASL: 110-240 V AC, 50/60 Hz, 2.0-1.0 A
Optics	

Description	Specification
Optical Platform	The system utilizes a fixed optical assembly configured with up to three spatially separated laser beams, with laser delays automatically adjusted during instrument Daily QC.
Optical Configurations	The following optical configurations are provided: - One-laser configuration: 488nm/50mW (Blue) - Two-laser configuration (B/R): 488nm/50mW (Blue), 640nm/80mW (Red) - Two-laser configuration (V/B): 405nm/100mW (Violet), 488nm/50mW (Blue) - Three-laser configuration: 405nm/100mW (Violet), 488nm/50mW (Blue), and 640nm/80mW (Red)
Beam Geometry	The system utilizes flat-top laser beam profile with narrow vertical beam height optimized for small particle detection.
Emission Collection	The system utilizes fused silica cuvette coupled to high NA lens for optimum collection efficiency to optical fibers.
Forward and Side Scatter Detection	The system utilizes a high-performance semiconductor detector with 488-nm bandpass filter for FSC detection, and two high-performance semiconductor detectors with 405-nm (only on systems containing Violet laser) and 488-nm bandpass filters for SSC detection.
Fluorescence Detection	The system utilizes a high sensitivity Coarse Wavelength Division Multiplexing (CWDM) semiconductor array per laser enabling more efficient spectrum capture in the 420-829 nm range. No filter changes are required for any fluorochrome excited by the on-board lasers.
Optical Detector Modules	The system utilizes the following optical detector modules: - Violet detector array (available in two- or three-laser configurations): 16 channels unevenly spaced bandwidth from 420-829 nm - Blue detector array: 14 channels unevenly spaced bandwidth from 498-829 nm - Red detector array (available in two- or three-laser configurations): 8 channels unevenly spaced bandwidth from 652-829 nm
Fluidics	
Fluidics Tasks	The vacuum driven fluidics module is capable of performing these fluidics tasks: Long Clean, SIT Flush, Purge Filter, Clean Flow Cell, Fluidics Shutdown.
Sample Flow Rates	The system is operable at the pre-set sample flow rates of: Low 15 µL/min, Medium 30 µL/min, High 60 µL/min
Fluidic Reservoirs	The system includes 4L sheath and waste tanks with built-in level-sensing and is compatible with 20L sheath and waste cubitainers.
Electronics	

Description	Specification
Signal Processing	The system is capable of digital signal processing with automatic window gate adjustment, 22-bit 6.5 log decades, threshold using any single parameter or combination of parameters.
Pulse Shape Parameters	The system provides Pulse Area and Height for every parameter, Width for scatter parameters and one fluorescence parameter for each laser.
Sample Input	
Compatible Tubes	The system is compatible with 12x75 mm polystyrene and polypropylene tubes.
Compatible ASL Carriers	The system is compatible with: U-, V-, and flat-bottom standard 96-well plate (polypropylene, polystyrene), 1-mL and 2-mL 96-deep well plate (polypropylene, polystyrene), and 40-tube rack
Minimum Volume	The minimum sample volume is 50 mL for all supported plates and tubes.
Transition Time	The switch from plate to tube rack acquisition takes less than 5 minutes.
Instrument Performance	
Fluorescence Linearity	The system fluorescence linearity defined as the correlation coefficient between known MESF and the measured fluorescence intensity is: $R \geq 0.98$ for FITC, PE, APC, and Pacific Blue channels
Forward and Side Scatter Resolution	The system has Forward and Side Scatter resolution optimized for resolving lymphocytes, monocytes, granulocytes and platelets from red blood cells.
Forward Scatter Resolution	The system is capable of resolving 1 μm microparticles from the debris by relative size.
Sample Carryover	The system sample carryover is: 1. < 0.5% in tube acquisition mode 2. < 0.3% in ASL Default acquisition mode for plate and tube rack 3. < 0.1% in ASL Low Carryover acquisition mode for plate and tube rack 4. < 1% in ASL High Throughput acquisition mode for plate and tube rack
Throughput Time ASL-Tube Rack	Acquisition of a 40-tube rack at standard 6C TBNK assay stopping rule of 3K lymphocytes with normal whole blood samples is completed: 1. < 20 minutes in High Throughput acquisition mode 2. < 27 minutes in Default acquisition mode
Throughput Time ASL-Plate	Acquisition of a 96-well standard or 96-deep well plate is completed: 1. < 27 minutes at 7 μL /well acquisition volume in High Throughput mode 2. < 22 minutes at 2 seconds/well acquisition time in High Throughput mode 3. < 48 minutes at 15 μL /well acquisition volume in Default mode

Description	Specification
ASL Dead Volume	ASL dead volume is: 1. < 10 mL for standard and deep well plates 2. < 40 mL for 12x75 mm tubes in 40-tube rack
ASL Sample Consumption Volume	ASL sample consumption volume is: 1. < 35 mL for 96-well flat-bottom plate 2. < 5 mL for all other supported plates and tubes
ASL Mixing	The Default mix settings is optimized for resuspending PBMC samples: 1. The %CV <10% for lymphocyte and granulocyte cell counts and percentage across all 96 wells within a plate. 2. The %CV <20% for Hela cell counts across all 96 wells within a plate.
ASL vs. Manual Tube Acquisition Equivalency	Between ASL and manual tube acquisition, for the lymphocyte subset percentage and cell counts in 6 Color TBNK assay samples: 1. The % difference of the mean between a 40-tube rack and manual tube is within $\pm 10\%$. 2. The % difference of the mean between a 96-well plate and manual tube is within $\pm 10\%$.
ASL Acquisition Precision	With ASL acquisition, for the lymphocyte subset percentage and cell counts in 6 Color TBNK assay samples: 1. The %CV from all 40 tubes in a tube rack is $\leq 10\%$. 2. The %CV from all 96 wells in a 96-well plate is $\leq 10\%$. 3. The difference between the 1st tube and last tube of a 40-tube rack is within $\pm 20\%$ of the 1st and last tube values respectively, calculated from the linear regression curve of all tubes. 4. The difference between the 1st well and last well of a standard 96-well plate is within $\pm 20\%$ of the 1st and last well values respectively, calculated from the linear regression curve of all wells.
Other Performance Indices	
Ploidy Analysis Linearity	For cell cycle analysis of diploid cells, the mean fluorescence intensity ratio of G2/M to G0G1 is between 1.95 and 2.05.
Surface Marker Detection Accuracy	For IVD control blood samples, the CD ₃ , CD ₄ , CD ₈ , CD ₁₆ /CD ₅₆ , and CD ₁₉ lymphocyte subset percentages are within the manufacturer provided ranges.
Surface Marker Detection Reproducibility	On repeated testing of control blood samples, the reproducibility for CD ₃ , CD ₄ , CD ₈ , CD ₁₆ /CD ₅₆ , CD ₁₉ lymphocyte subset percentages meets: 1. The CV $\leq 8\%$ when the subset percentage is greater than or equal to 30% or 2. The CV $\leq 15\%$ when the subset percentage is less than 30%

Description	Specification
Instrument Stability	When ambient temperature change does not exceed 5% of the set temperature, the fluctuation range within an 8-hour period is no greater than $\pm 10\%$ for detected forward scatter and peak fluorescence channels among all fluorescence channels.

Supplies and Replacement Parts

Item	Part Number	Description
4-L Tank	N4-00124	Fluidics bottle for sheath and waste
4L Supply Cap Assy	N7-32015	Lid for 4L sheath tank with tubing
4L Waste cap Assy	N7-32014	Lid for 4L waste tank with tubing and liquid level sensor
20L Cubitainer Adapter Assy	N7-32016	Tubing and caps to fit 20L waste and sheath containers
Reservoir Holder Assy	N7-32019	Holder for 4L sheath and waste tanks
SpectroFlo® QC Beads 2000 Series, 2mL	B7-20002	Daily QC beads to validate instrument performance
Sheath filter assembly (0.2 µm)	N7-22006	0.2-µm sheath filter assembly with quick-connect fittings and manifold
Sheath filter assembly (0.4 µm)	N7-22016	0.4-µm sheath filter assembly with quick-connect fittings and manifold
Sample line	N7-22014	SIT tubing assembly
Sheath filter bypass (Long Clean tubing)	N7-22010	Replaces the sheath filter during a Long Clean
25-mm disk filter (20 µm)	N4-00047	Filter at sheath tubing intake for cubitainer and 4-L sheath tank cap
Cytek Sheath Fluid	TNB-4600-L600	Cytek Biosciences Sheath Fluid (30X)
Cytek Shutdown Fluid	TNB-4601-L1500	Shutdown Fluid is a ready-to-use 1X solution for long-term storage of Cytek Bioscience analyzers

A

accumulator vessel 18, 19

acquisition

controls 49

experiments 48

preferences 75

adding

fluorescent tags to library 68

labels 69

new user account 86

alarm ranges 45

Analyzing Data Offline 101

area scaling factor 51

B

benchmark reference controls 43

C

cleaning

external surfaces 139

flow cell 137

SIT 136

cytometer

cleaning surfaces 139

fluid containers 18, 29–30

fluidics system components 18–19

lasers 20

optics 19

overview 16

power button 16

sheath filter 19

shutdown 32

SIT/SIP 17

D

daily QC

performing 33–36

report 36

data

formats 47

raw 47, 101

unmixed 47, 101, 102, 105

degasser 18

deleting

user account 87

E

editing

fluorescent tag properties 68

gate properties 54

plot properties 53

user account 87

experiment

about 23

creating new 58–64

default 23, 27

display 48–55

exporting 24

setup overview 48

templates 24, 73

toolbar 52

export

experiment 24

QC report 36

F

FCS files

formats 24, 47

raw data 101

storage preferences 82

unmixed data 101, 102, 105

flow cell, cleaning 137

flow meter 141

fluid containers 18, 29–30

fluidics lines

inspecting 139

fluidics system

components 18–19

decontaminating 137

fluorescent tag

adding new 68

editing properties 68

groups 38

library 67

live unmixing 97

running reference controls 38

- selecting in experiment 58
- font preferences 81

G

- gain

- adjusting 51
 - See also* user settings
 - settings 45

- gates

- adjusting, post-acquisition unmixing 101, 105
 - editing properties 54
 - preferences 79

K

- keywords 69

L

- labels

- adding 69
 - library 68

- lasers

- detectors 20
 - scaling factor 51

- Levey-Jennings

- tab for performance tracking 33
 - tracking 44

- library

- experiment templates 73
 - fluorescent tags 67
 - keywords 69
 - labels 68
 - loader settings 72
 - QC beads 67
 - user settings 71
 - worksheet templates 72

- loader

- calibrating the loader Stage and plate 143
 - connections 145
 - daily workflow 127
 - lifting and carrying instructions 145
 - maintenance 143
 - set up experiment 127
 - settings 72
 - settings and controls 125
 - shut down 132
 - startup 127
 - troubleshooting 150
 - uninstalling power 145

- long clean 137

M

- maintenance

- scheduled 135
 - unscheduled 135

- modules, software 22

N

- notification preferences 82

O

- optics 19

- compartment 16

- overview

- cytometer 16
 - Northern Lights system 15
 - software 22
 - unmixing 94

P

- password, resetting 87

- plenum 18

- plots

- preferences 78
 - types and properties 52

- post-acquisition unmixing

- in acquisition module 101–102
 - in analysis module 102, 103–105

- power button 16

- preferences

- acquisition 75
 - fonts 81
 - gates 79
 - notifications 82
 - plot 78
 - QC setup 84
 - statistics 80
 - storage 82
 - worksheet 77

Q

- QC beads

- library 67
 - running 33–36

- QC reports

- exporting 36
 - pass/fail criteria 36

- QC setup

- preferences 84

R

- reference controls 38
 - comparing to benchmarks 98
 - running 38–42
 - selecting benchmarks 43
 - updating 43

- replacing
 - sheath filter 139
 - SIT 140–143

S

- safety
 - biological 12
 - electrical 12
 - general information 11
 - symbols 9
- sample injection port (SIP) 17
- sample injection tube (SIT) 17
 - compartment 16
 - flushing 136
 - flushing, loader 127
 - replacing 140–143
- settings
 - user 45
- sheath container 18
 - filling 29–??
- sheath filter 19
 - purging 136
 - replacing 139
- sheath plenum 18
- shutdown procedure 32
- shutting down system 32
- signal measurements 51
- SIP
 - See sample injection port
- SIT
 - See sample injection tube
- software
 - FCS file formats 47
 - modules 22
 - overview 22
- spillover 93
 - adjusting 106
- statistics
 - creating 54
 - preferences 80
- system
 - overview 15
 - shutdown 32
 - starting up 31

T

- technical support information 13
- threshold 51
- toolbar, experiment 52
- tracking
 - Levey-Jennings 44
- troubleshooting 147–150

U

- unmixing 93
 - label/lot specific 76
 - live 96–101
 - overview 94
 - post-acquisition 101–105
- user settings 45
 - library 71
- users
 - adding new account 86
 - editing account 87
 - removing account 87
 - resetting user password 87
 - use time on system 88

V

- virtual filters 106–110

W

- waste container 18
 - emptying 30
- worklist
 - add tasks 112
 - arrange sample locations in a carrier 116
 - create or open 111
 - data collection 116
 - export report 117
 - lot specific reference controls 115
 - reference controls 113
 - unmix 114
- worksheet
 - about 26
 - preferences 77
 - selecting 52
 - templates 72



 Cytek Biosciences, Inc
47215 Lakeview Blvd.
Fremont, CA 94538, USA
Phone: 1.877.922.9835 (option 1)

Cytek Biosciences B.V.
Paasheuvelweg 25, Tower C5
1105 BP Amsterdam, The Netherlands
Phone: +31 (0) 20 765 3440

cytekbio.com



Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands